BIOSYNEX Bladder Tumor Antigen

REF 1110004

FRANÇAIS (FRENCH)3

ENGLISH.....7

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MANUFACTURER INFORMATION

BIOSYNEX S.A.

22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN – France
Standard : +33 3 88 78 78 87
www.biosynex.com

Contact France : +33 3 88 77 57 00

service.clients@biosynex.com

Contact other countries : +33 3 88 77 57 52

sales@biosynex.com

After-sales service: +33 3 88 77 57 25


tech.support@biosynex.com

CH	REP
----	-----

BIOSYNEX SWISS S.A.

Route de Rossemaison 100
2800 DELEMONT - Switzerland

NON-STANDARD SYMBOLS

	FR	EN	DE
	Diluant	buffer	Verdünnungspuffer

BIOSYNEX Bladder Tumor Antigen

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13. BIBLIOGRAPHIE

1. Badrinath R. Konety and Robert H. Getzenberg: Urine based markers of urological malignancy. The Journal of Urology, Vol. 165, 600-611, Feb, 2001
2. M-P Raitanen, E Kaasinen, E Rintala, E Hansson, P Nieminen, R Aine, TLJ Tammela1 and The Finn Bladder Group*: Prognostic utility of human complement factor H related protein test. British Journal of Cancer (2001) 85(4), 552-556
3. Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. Symposium on Diagnostic Accuracy of Cytologic Technics 8:186-193, 1964.
4. Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. Cancer 60:1423-1427, 1987.
5. Robert-Koch Institut (Stand 1/2009) please refer to: http://www.rki.de/clin_100/nn_203956/DE/Content/GBE/Dachdok_Krebs/Broschuere/Lokalisationen/C67_D090_D414,templateId=raw,property=publicationFile.pdf/C67_D090_D414.pdf
6. Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA stat Test. Urology 50(3):349 - 353, 1997.
7. Raitanen M.-P., Marttila T., et al: The Bard BTA stat Test in Monitoring of Bladder Cancer. The Journal of Urology 157: 28, 1997.
8. Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Carcinoma of the Bladder. Clinical Cancer Research 4:2511-2520, 1998.
9. Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. Proceedings of the American Association for Cancer Research 39: 263, 1998.

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BIO SYNEX Bladder Tumor Antigen

A RAPID TEST FOR THE QUALITATIVE DETECTION OF BLADDER TUMOR ANTIGEN IN URINE.

In vitro diagnostic medical device for professional use.

REF 1110004

1. INTENDED USE

The BIOSYNEX Bladder Tumor Antigen test is a simple one step lateral flow immunoassay for the qualitative detection of Bladder Tumor Marker in urine.

The BIOSYNEX Bladder Tumor Antigen test can be used as an aid to initiate or attend therapeutic treatments by physicians.

2. CLINICAL SUMMARY

2-3% of all human malign tumors are bladder cancer tumors. After prostate carcinoma bladder cancer tumor comes in second of all tumors of the genitourinary system. Younger people are not so heavy concerned, the risk is increasing with age. Men are having a three times higher risk for this disease as women.

Certain exposures such as smoking, some chemicals and drugs stimulate the development of bladder cancer.

Established diagnostic methods to diagnose bladder tumors are urine cytology and cystoscopy. The disadvantage of urine cytology is a variable sensitivity; especially tumors in their early stages are often missed. Cystoscopy is not without risk because of its invasive nature furthermore will it only identify visible tumors.

Therefore, a marker for diagnosis of bladder cancer has been searched in recent time. Bladder tumor protein concentration is elevated in urine samples of patients with bladder cancer. Elevated level of this protein is detected by the BIOSYNEX Bladder Tumor Antigen test.

Several studies have confirmed that bladder tumor protein is a useful marker for the detection of bladder cancer. BIOSYNEX Bladder Tumor Antigen is an easy-to-use and patient friendly rapid test to diagnose and monitor bladder cancer.

3. TEST PRINCIPLE

The BIOSYNEX Bladder Tumor Antigen test employs a unique combination of monoclonal and polyclonal antibodies to selectively identify Bladder Tumor Protein in test samples with a high degree of sensitivity.

Bladder Tumor Protein-Antibodies are immobilized on the test region of the membrane. During the test, the specimen reacts with other Bladder Tumor Protein-antibodies conjugated to colored particles and precoated onto the sample pad of the test. The mixture then migrates through the membrane by capillary action and interacts with reagents on the membrane. If there is sufficient Bladder Tumor Protein in the specimen, a colored line will form at the test region of the membrane. The presence of this colored line indicates a positive result, while its absence indicates a negative result. The appearance of a colored line at the control region serves as a procedural control, indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. REQUIRED MATERIAL

Material provided:

- Cassettes and Pipettes packed in individual pouch containing desiccant bag

Material required but not provided:

- Timer
- Specimen collection container for urine samples collection

5. PRECAUTIONS

- For professional in vitro diagnostic use only.
- Do not use after the expiration date.
- For optimal results, follow the procedure and storage conditions carefully.
- Samples and used materials must be considered as potentially infectious and treated as such. Disposed of used components according to the procedure and local regulations applicable to potentially infectious waste. Clean surface and any spilled liquid with an appropriate disinfectant.
- Do not use test if pouch is damaged.
- Do not use more than the required amount of liquid.
- Use only urine as liquid and no other one instead.
- Avoid cross-contamination of urine samples by using a new specimen collection container and specimen pipette for each urine sample.

- Follow Good Laboratory Practice when collecting the specimen and performing the test.
- Store and transport the test device always at 2-30°C (36°-86°F).
- Inadequate humidity and temperature conditions may adversely affect results.

6. KIT STORAGE

The kit should be stored at 2-30°C until the expiry date printed on the box and test must remain in the sealed pouch until use. **Do not freeze.**

Care should be taken to protect kit's components from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment or containers can lead to false results.


7. SAMPLE/CONTROLS COLLECTION AND STORAGE

- Specimen collection should not be performed during or within three days of a menstrual period for female.
- Only urine specimens should be used in this assay. Specimen should avoid contamination of toilet water.
- Specimens collected, if not tested the same day, can be refrigerated and stored at 4 to 8 °C for testing up to 3 days after collection.

8. TEST PROCEDURE

1. Bring tests and specimens to room temperature (15-30°C) before use.
2. Remove the test cassette from the foil pouch, and place it on a flat, dry surface. Label the test with patient or control identification. For best results, the assay should be performed within one hour.
3. Using the provided disposable pipette, transfer 2-3 drops of specimen to the specimen well (S) of the device (Figure below) and start the timer. Avoid trapping air bubbles in the specimen well (S), and do not add any solution to the result area.
4. As the test begins to work, color will migrate across the membrane. Wait for the colored line(s) to appear. The result should be read at 9 minutes.

Do not interpret test results after 10 minutes.

9. INTERPRETATION OF RESULTS

- POSITIVE*:**
- Presence of 2 distinct colored lines:** A control line appears at the level of the C zone and one colored line (even of weak intensity) appears at the level of the test T zone. A positive result indicates that Bladder Tumor Antigen has been detected.
- NEGATIVE:**
- Only one colored line appears at the level of the control line (C).**
- No line appears at the level of the test T zone. A negative result indicates that the specimen contains no Bladder Tumor Antigen
- INVALID:**
- No visible colored band at the level of control line C (whatever test line T apparition).**
- Results from a test with no control line must be discarded. Review the procedure and repeat it with a new test device. If the problem persists, contact your local distributor.

Notes:

- The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any color in the test region (T) should be considered positive. Note that this is a qualitative test, and cannot determine the concentrations of analytes in specimens.
- Once a positive result has been established (at 10 minutes), the result will not change. However, in order to prevent any incorrect results, the test result should not be interpreted after 10 minutes
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most frequent reasons for control line failure.

10. QUALITY CONTROL

- Internal procedural controls are included in the test. A colored line appearing in the control zone (C) ensures that sufficient specimen volume has been loaded and that the correct procedure has been followed by the operator.
- Good laboratory practices recommend the use of control materials to ensure proper kit performance. Control samples specific for this product are available separately.

11. PERFORMANCE

A. Clinical study

A comparison study has been performed with cytology. A total of 225 patients with biopsy confirmed bladder cancer took part in this study. Of the 225 patients, stages were determined on 223, and grades are determined on 210 patients.

Stage	N	Sensitivity
Ta	14	57%
T1	69	82.6%
T2,T3	140	86.4%
Total	223	83.4%
Grade	N	
1	47	38.2%
2	61	63.9%
3	102	89.2%
Total	210	70.5%

Clinical Sensitivity:

The overall sensitivity varies between 70.5 % to 83.4% depending on the stages and grades of the bladder tumor. The sensitivity of the test varies from 38.2% (Grade 1) to 89.2% (Grade 3).

N	True Negative
Bladder Positive	4
Bladder Negative	34
Total	38

The test has an overall specificity of 89.5% (34/38).

B. Analytical Specificity

Cross reactivity:

Positive and negative samples have been spiked with following substances:

<i>Acinetobacter calcoaceticus</i>	<i>Neisseria catarrhalis</i>
<i>Proteus vulgaris</i>	<i>Gardnerella vaginalis</i>
<i>Salmonella typhi</i>	<i>Neisseria meningitidis</i>
<i>Acinetobacter spp</i>	<i>Streptococcus faecalis</i>
<i>Staphylococcus aureus</i>	<i>Neisseria lactamica</i>
<i>Candida albicans</i>	<i>Streptococcus faecium</i>
<i>Neisseria gonorrhoeae</i>	<i>Pseudomonas aeruginosa</i>
<i>Escherichia coli</i>	<i>Trichomonas vaginalis</i>

None of the tested substances showed any cross-reaction with neither a clinically defined negative nor a positive specimen. Negative specimen samples with supplementation of potentially cross reacting substances gave consistently negative test results, whereas specimen samples positive to Bladder Tumor scored consistently positive.

Interferences:

Possible interfering substances were added to negative and positive urine samples with following amounts :

Acétamiophene 20 mg/dL	Ketones 40 mg/dL
Acetyl salicylic Acid 20 mg/dL	Mestranol 3 mg/dL
Ascorbic Acid 20 mg/dL	Nitrite 20 mg/dL
Atropine 20 mg/dL	Penicillin 40,000 U/dL
Bilirubin 60 mg/dL	Sodium Heparin 3 mg/dL
Caffeine 20 mg/dL	Lithium Heparin 3 mg/dL
Creatinine 20 mg/dL	LDL 200 mg/dL
Gentesic Acid 20 mg/dL	Triglycerides 150 mg/dL
Glucose 2000 mg/dL	

None of the tested substances showed any interference with neither a clinically defined negative nor a positive specimen. Negative specimen samples with supplementation of potentially interfering substances gave consistently negative test results, whereas specimen samples positive to Bladder Tumor scored consistently positive.

12. LIMITATIONS

- If a positive result is obtained with a patient taking aspirin, the test should be repeated after aspirin has been stopped for 7 days.
- A positive Bladder Tumor Antigen result indicate a likelihood of Bladder Cancer, however requires further confirmation.
- Negative results do not exclude bladder cancer. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- The assay is designed for use with human urine only. Due to too low ionic strength (e.g. absence of ions) of some samples, it might be possible, that false results occur.

13. LITERATURE

1. Badrinath R. Konety and Robert H. Getzenberg: Urine based markers of urological malignancy. The Journal of Urology, Vol. 165, 600-611, Feb, 2001
2. M-P Raitanen, E Kaasinen, E Rintala, E Hansson, P Nieminen, R Aine, TLJ Tammela¹ and The Finn Bladder Group*: Prognostic utility of human complement factor H related protein test. British Journal of Cancer (2001) 85(4), 552–556
3. Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. Symposium on Diagnostic Accuracy of Cytologic Technics 8:186-193, 1964.
4. Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. Cancer 60:1423-1427, 1987.
5. Robert-Koch Institut (Stand 1/2009) please refer to: http://www.rki.de/clin_100/nn_203956/DE/Content/GBE/Dachdok_Krebs/Broschuere/Lokalisationen/C67_D090_D414,templateId=raw,property=publicationFile.pdf/C67_D090_D414.pdf
6. Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA stat Test. Urology 50(3):349 - 353, 1997.
7. Raitanen M.-P., Marttila T., et al: The Bard BTA stat Test in Monitoring of Bladder Cancer. The Journal of Urology 157: 28, 1997.
8. Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Carcinoma of the Bladder. Clinical Cancer Research 4:2511-2520, 1998.
9. Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. Proceedings of the American Association for Cancer Research 39: 263, 1998.

14. MODIFICATION HISTORY

IFU version	Amended paragraph	Details of the changes
V1	All	Multilingual format and transition to V1 Update to biosynex format
	Content	Addition of the table of contents
	14 : History	Addition of this paragraph
V2	4 : Required material	Removal of instruction for use from the kit

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13. LITTERATUR

1. Badrinath R. Konety and Robert H. Getzenberg: Urine based markers of urological malignancy. The Journal of Urology, Vol. 165, 600-611, Feb, 2001
2. M-P Raitanen, E Kaasinen, E Rintala, E Hansson, P Nieminen, R Aine, TLJ Tammela¹ and The Finn Bladder Group*: Prognostic utility of human complement factor H related protein test. British Journal of Cancer (2001) 85(4), 552-556
3. Umiker W.: Accuracy of Cytologic Diagnosis of Cancer of the Urinary Tract. Symposium on Diagnostic Accuracy of Cytologic Technics 8:186-193, 1964.
4. Badalament R.A., Hermansen D.K., Kimmel M., Gay H., Herr H.W., Fair W.R., Whitmore W.F., Jr., Melamed M.R.: The Sensitivity of Bladder Wash Flow Cytometry, Bladder Wash Cytology, and Voided Cytology in the Detection of Bladder Carcinoma. Cancer 60:1423-1427, 1987.
5. Robert-Koch Institut (Stand 1/2009) please refer to: [http://www.rki.de/cin_100/nn_203956/DE/Content/GBE/Dachdok Krebs/Broschuere/Lokalisationen/C67_D090_D414,templateId=raw,property=publicationFile.pdf/C67_D090_D414.pdf](http://www.rki.de/cin_100/nn_203956/DE/Content/GBE/Dachdok/Krebs/Broschuere/Lokalisationen/C67_D090_D414,templateId=raw,property=publicationFile.pdf/C67_D090_D414.pdf)
6. Sarosdy M.F., Hudson M.A., et al: Improved Detection of Recurrent Bladder Cancer Using the Bard BTA stat Test. Urology 50(3):349 - 353, 1997.
7. Raitanen M.-P., Marttila T., et al: The Bard BTA stat Test in Monitoring of Bladder Cancer. The Journal of Urology 157: 28, 1997.
8. Kinders R., Jones T., et al: Complement Factor H or a Related Protein Is a Marker for Transitional Cell Carcinoma of the Bladder. Clinical Cancer Research 4:2511-2520, 1998.
9. Corey M., Kinders R., et al: Factor H Related Proteins Are Upregulated In Bladder Cancer. Proceedings of the American Association for Cancer Research 39: 263, 1998

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