BIOSYNEX Strep A

REF 1030005

FR (FRANÇAIS / FRENCH)	1
EN (ENGLISH)	6
DE (DEUTSCH / GERMAN)	10
ES (ESPAÑOL / SPANISH)	15
IT (ITALIANO / ITALIAN)	20
PT (PORTUGUES / PORTUGUESE)	25
CS (ČESKY / CZECH)	29
NO (NORSK / NORWEGIAN)	33
DA (DANSK / DANISH)	37
FI (SUOMALAINEN / FINNISH)	41
SV (SVENSKA / SWEDISH)	45

CONTENTS

MANUFACTURER INFORMATION

BIOSYNEX SWISS S.A.
Route de Rossemaison 100
2800 DELEMONT - Switzerland
Standard : **+33 3 88 78 78 87**
www.biosynex.com

Contact France : **+33 3 88 77 57 00**
service.clients@biosynex.com
Contact other countries : **+33 3 88 77 57 52**
sales@biosynex.com
After-sales service: **+33 3 88 77 57 25**
tech.support@biosynex.com

EU	REP
----	-----

BIOSYNEX S.A.
22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN - France

NON-STANDARD SYMBOLS

	FR	EN	DE	ES	IT	PT	CS	NO	DA	FI	SV
DIL	Diluant	Buffer	Verdünnungspuffer	Tampón	Tampone di diluizione	Tampão	Pufr	Buffer	Buffer	Puskuri	Buffertlösning

CONTENTS

BIOSYNEX Strep A

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BIOSYNEX Strep A

RAPID TEST FOR THE QUALITATIVE DETECTION OF GROUP A STREPTOCOCCI FROM OROPHARYNGEAL SWABS.

In vitro diagnostic medical device for professional use.

REF 1030005

1. INTENDED PURPOSE

BIOSYNEX Strep A is a rapid lateral flow immunoassay for the qualitative, presumptive detection of Group A Streptococcus antigens from oropharyngeal specimens. This kit is intended for use as an aid in the diagnosis of Strep A infections.

2. CLINICAL SUMMARY

Streptococcus pyogenes are non-motile gram-positive cocci, which can colonize different parts of the human body and cause serious infections.

Beta-haemolytic Group A Streptococci (*Streptococcus pyogenes*) are the main cause for infections of the upper respiratory tracts like tonsillitis, pharyngitis, and other respiratory infections, moreover impetigo, endocarditis, puerperal sepsis, meningitis, and arthritis. Left untreated, these infections can lead to serious complications, including rheumatic fever and peritonsillar abscesses. An early diagnosis and treatment of Group A Streptococcal pharyngitis has been shown to reduce the severity of symptoms and further complications, such as rheumatic fever and glomerulonephritis.

Conventional methods for detecting Strep A infections are dependent on isolation and subsequent identification of the organism, and often require 24-48 hours.

BIOSYNEX Strep A is a rapid test for the detection of Strep A antigens in throat providing results within 5 minutes. It helps the medical practitioner make a rapid diagnosis and set up an immediate and selective therapy. BIOSYNEX Strep A utilizes antibodies specific for whole cell Lancefield Group A Streptococcus for the sensitive detection of Strep A antigens in oropharyngeal specimens.

3. TEST PRINCIPLE

BIOSYNEX Strep A is a qualitative, lateral flow immunoassay for the detection of Strep A carbohydrate antigen in oropharyngeal specimens.

In this test, antibodies specific to Strep A antigens are immobilized in the test line area (T). During the test, the antigens extracted from the swab specimen are captured by Strep A-specific antibodies, which are adhered to colored particles. The mixture migrates along the membrane and the antigen-antibody-particle complex binds to the specific antibody in the test line area (T). The agglomeration of complexes creates a colour line in the test line area (T).

The appearance of the colour line in the test (T) line area indicates a positive result, while its absence indicates a negative result. A red line should always appear in the control line area (C). It serves as a procedural control and indicates that the test ran correctly and that the result is valid.

4. REAGENTS

The test devices include Strep A antibody coated colored particles and Strep A antibodies coated on the membrane.

5. REQUIRED MATERIAL

Material provided:

- Test cassettes
- Extraction solution 1 (2 M NaNO₂)

Danger

Harmful if swallowed

- Extraction solution 2 (Citric acid 0.027 M)
- Positive control (non-viable Strep A; 0.09 % NaN₃)
- Negative control (non-viable Strep C; 0.09 % NaN₃)
- Extraction tube with dropper caps
- Sterile swabs
- Reagent holder

Material required but not provided:

- Timer

6. PRECAUTIONS

- For single use only.
- Do not use components after stated expiration date (see pouch and box label).
- Do not use test or swab, if pouch is damaged.

- Handle all samples as infectious products. Observe established precautions against microbiological hazards during the procedure and follow standard procedures for sample disposal.
- Samples and used materials must be considered as potentially infectious and treated as such. Dispose of used components according to the procedure and local regulations applicable to potentially infectious waste. Clean surface and any spilled liquid with an appropriate disinfectant.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Humidity and high temperature can adversely affect results.
- Extraction solution 1: Wash hands thoroughly after handling. Do not eat, drink or smoke when using this product. IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell. Rinse mouth. Dispose of contents/container in accordance with local/regional/national/international regulations.
- Do not use more than the required amount of liquid.
- Do not spill the specimens into the reaction area. Do not touch the reaction area of the device to avoid contamination.
- Do not mix reagents from different lots. Do not mix reagent bottle caps.
- Avoid cross-contamination of specimens by using a new extraction tube and specimen pipette for each specimen.
- Use only Dacron or Rayon tipped sterile swabs with plastic shafts such as those provided. Do not use calcium alginate, cotton tipped, or wooden shafted swabs.
- Positive and negative controls contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these solutions always flush with copious amounts of water to prevent azide buildup.
- The potentially infectious materials (e. g. antibodies) or other components of the test (chemicals) do not constitute any danger if test is used according to instructions.
- Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.
- Safety datasheet available upon request.

7. KIT STORAGE

The kit should be stored at 2-30°C. The test is stable through the expiry date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE.

8. SPECIMEN COLLECTION AND STORAGE

1. Using the sterile swab provided in the test kit, collect the specimen on the posterior pharynx, tonsils and other inflamed areas. Avoid touching the tongue, cheeks and teeth with the swab.
2. Testing should be performed immediately after collecting the sample. Swab specimens can be stored in a clean, dry plastic tube for up to 8 hours at room temperature (15-30°C) or 72 hours at 2-8°C.
3. If a confirmatory culture is preferred, lightly roll the swab on a Group A (GAS) selective blood agar plate before using the swab provided with BIOSYNEX Strep A.

9. TEST PROCEDURE

Bring tests, Extraction solution s, swab specimens, and/or external controls to room temperature (15-30°C) before testing.

1. Remove the test from its sealed pouch and use it within the hour. For best results, the assay should be performed within immediately after opening the sealed pouch.
2. Hold the Extraction solution 1 bottle vertically and add **4 full drops** (approximately 240 µL) to an extraction tube. Extraction solution 1 is red in colour.
3. Hold the Extraction solution 2 bottle vertically and add **4 full drops** (approximately 160 µL) to the extraction tube. Extraction solution 2 is colourless. Mix the solution by carefully swirling the extraction tube. The addition of Extraction solution 2 to Extraction solution 1 changes the colour of the solution from red to yellow.
4. Add immediately the swab into the extraction tube; agitate the swab vigorously 15 times while pressing the head against the bottom of the tube to release the antigen in the swab.
5. Press the swab against the wall of the tube and squeeze the bottom of the tube while removing the swab. Make sure to have as much solution as possible remain in the tube. Discard the swab and let incubate the solution for at least 1 minute. The test sensitivity increases with prolonged time of extraction. Therefore, it is recommended to extract for at least 5 min, because it will affect the test result at concentrations near the detection limit.
6. Place the dropper cap on the extraction tube. Place the test device on a clean and level surface. Add **3 drops** of the solution (approximately 200 µL) from the extraction tube to the sample well. Avoid trapping air bubbles in the specimen well (S) and do not add any liquid to the reaction area. Start the timer as the test starts to run.
7. **Interpret results after 5 minutes.** Do not interpret any results after more than 10 minutes.

10. RESULTS INTERPRETATION

POSITIVE:


2 lines appear. One line appears in the control line area (C) and one line in the test line area (T). A positive result indicates that Strep A has been detected.

NOTE: The intensity of colour in the test area (T) may vary depending on the concentration of Strep A present in the specimen. Therefore, any shade of colour in the test area (T) should be considered positive.

NEGATIVE:

One line appears in the control line area (C). No line appears in the test line area (T). A negative result indicates that no Strep A antigen is present in the specimen or that it is below the detection level of the test device. The specimen should be cultured to confirm the absence of a Strep A infection. If clinical symptoms are not consistent with the test results, obtain another specimen for culture.

INVALID:

Control line fails to appear. Insufficient specimen volume, expired test components or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

11. QUALITY CONTROL

Internal controls:

A procedural internal control is included in the test. A coloured line appearing at the control line (C) is considered a positive procedural internal control, confirming sufficient sample volume and correct procedure.

External controls:

It is recommended to perform a positive and negative external control for every kit, and as deemed necessary by internal laboratory procedures. External positive and negative controls are supplied with the kit. Alternatively, other Group A and non-Group A Streptococcus reference strains may be used as external positive controls. Please note, that commercial controls may contain interfering preservatives; therefore, other commercial controls are not recommended.

Procedure for External Quality Control Testing:

1. Add **4 full drops** of Extraction solution 1 and **4 full drops** of Extraction solution 2 into an extraction tube.
2. Add **1 full drop** of positive or negative control solution into the tube, holding the bottle vertically upside down.
3. Place a clean swab into the extraction tube and mix the solution by rotating the swab at least 15 times. Leave the swab in the extraction tube for 1 minute. Then express the liquid from the swab head by rolling the swab against the inside of the extraction tube and squeezing the extraction tube as the swab is withdrawn. Discard the swab.
4. Continue with Step 6 of the „TEST PROCEDURE“ section. If the control does not yield the expected results, do not use the test results. Repeat the procedure for external quality control testing or contact your distributor.

12. PERFORMANCE

Diagnostic Sensitivity and Specificity

Using three medical centres for evaluation, a total of 525 samples were collected from patients showing symptoms of pharyngitis. Each swab was rolled onto a sheep blood agar plate, and then tested with BIOSYNEX Strep A.

		Culture		Total
		Positive	Negative	
BIOSYNEX Strep A	Positive	117	11	128
	Negative	7	390	397
	Total	124	401	525

Relative sensitivity 94.4 % [95%CI*: 88.7%-97.7%]

Relative specificity 97.3 % [95%CI*: 95.1%-98.6%]

Correlation 96.6 % [95%CI*: 94.6%-98.0%]

*Confidence interval

If results were aligned to frequency of Strep A in culture the following results were obtained.

Positive Culture Classification	Strep A Rapid Test/Culture	% Compliance
Rare	10/12	83,3%
1+	20/22	90,9%
2+	18/20	90,0%
3+	31/32	96,9%
4+	38/38	100,0%

Analytical Sensitivity

The analytical sensitivity of the assay is 1×10^5 bacteria / swab. 8 different strains of Strep A were tested and all showed weak positive results at this concentration.

High Dose Hook

No High Dose Hooke Effect was observed up to a concentration of 1.0×10^{10} bacteria / swab. This indicates that the measurement range is at least 1.0×10^5 to 1.0×10^{10} bacteria / swab.

Inter-/Intra LOT variances

Three independent LOT were tested in 5 fold determinations with negative controls and low positive controls. All assays showed the expected results independent of LOT and determination. From this it can be concluded that Inter-/Intra LOT variance is low.

Cross Reactivity

The following organisms were tested at 1.0×10^7 organisms per test and were all found to be negative when tested with BIOSYNEX Strep A. No mucoid-producing strains were tested.

Group B Streptococcus	Staphylococcus aureus	Neisseria gonorrhea
Group C Streptococcus	Candida albicans	Neisseria meningitides
Group F Streptococcus	Corynebacterium diphtheria	Neisseria sicca
Group G Streptococcus	Branhamella catarrhalis	Neisseria subflava
Streptococcus pneumoniae	Serratia marcescens	Pseudomonas aeruginosa
Streptococcus sanguis	Klebsiella pneumoniae	Enterococcus faecalis
Streptococcus mutans	Bordetella pertussis	
Staphylococcus epidermidis	Hemophilus influenza	

13. EXPECTED VALUES

Approximately 15% of pharyngitis in children ages 3 months to 5 years is caused by Group A beta haemolytic Streptococcus. In school-aged children and adults, the incidence of Strep throat infection is about 40%. This disease usually occurs in winter and early spring time in temperate climates.

14. LIMITATIONS

- The test should be used for the detection of Strep A antigen in oropharyngeal swab samples. The Strep A antigen concentration cannot be determined by this test.
- This test indicates the presence of Strep A antigen in the sample from both viable and nonviable Group A Streptococcus bacteria.
- The test does not differentiate asymptomatic carriers of Group A Streptococcus from those with a symptomatic infection. A negative result should be confirmed by culture. A negative result may be obtained if the concentration of the Strep A antigen present in the swab specimen is not adequate or is below the detectable level of the test.
- Excess blood or mucus on the swab specimen may interfere with test performance and may yield a false result. When collecting the pharyngeal swab specimen avoid contact with tongue, cheeks, teeth or any bleeding areas of the mouth.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single rapid test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- The accuracy of the test depends on the quality of the swab specimen. False negatives may result from improper specimen collection or storage. A negative result may also be obtained from patients at the onset of the disease due to low antigen concentration.
- Respiratory infections, including pharyngitis, can be caused by streptococci from serotypes other than Group A, as well as other pathogens.

15. LITERATURE

1. Murray, P.R., et al. Manual of Clinical Microbiology, 6th Edition, ASM Press, Washington D.C., 1995, p. 299-307.
2. Webb, KH. Does Culture Confirmation of High-sensitivity Rapid Streptococcal Tests Make Sense? A Medical Decision Analysis. Pediatrics (Feb 1998), 101:2, 2.
3. Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. Diagnosis and Management of Group A Streptococcal Pharyngitis. Clinical Infectious Diseases (1997), 25: 574-83.
4. Needham CA, McPherson KA, Webb KH. Streptococcal Pharyngitis: Impact of a High-sensitivity Antigen Test on Physician Outcome. Journal of Clinical Microbiology (Dec 1998), 36: 3468-3473.
5. Shea, Y.R., Specimen Collection and Transport, Clinical Microbiology Procedures Handbook, Isenberg, H.D., American Society of Microbiology, Washington D.C., 1.1.1- 1.1.30, 1992.
6. Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. Group A beta-hemolytic streptococcal pharyngitis in preschool children aged 3 months to 5 years. Clinical Pediatrics (June 1999), 38: 357-360.
7. Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA. Group A Streptococcal Pharyngitis in Adults 30 to 65 years of age. Southern Medical Journal (May 1999), 491- 492.

16. MODIFICATION HISTORY

IFU version	Amended paragraph	Details of the changes
V1	All	Multilingual format and transition to V1
	Content	Addition of the table of contents
	16 : Modification history	Addition of this paragraph
V2	Content	Changed the EC Rep symbol to the EU Rep symbol
	5: Required material	- Removed the CE number of the sterile swabs - Removed the package insert from the material provided

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15. LITTERATUR

1. Murray, P.R., et al. Manual of Clinical Microbiology, 6th Edition, ASM Press, Washington D.C., 1995, p. 299-307.
2. Webb, KH. Does Culture Confirmation of High-sensitivity Rapid Streptococcal Tests Make Sense? A Medical Decision Analysis. Pediatrics (Feb 1998), 101:2, 2.
3. Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. Diagnosis and Management of Group A Streptococcal Pharyngitis. Clinical Infectious Diseases (1997), 25: 574-83.
4. Needham CA, McPherson KA, Webb KH. Streptococcal Pharyngitis: Impact of a High-sensitivity Antigen Test on Physician Outcome. Journal of Clinical Microbiology (Dec 1998), 36: 3468-3473.
5. Shea, Y.R., Specimen Collection and Transport, Clinical Microbiology Procedures Handbook, Isenberg, H.D., American Society of Microbiology, Washington D.C., 1.1.1- 1.1.30, 1992.

6. Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. Group A beta-hemolytic streptococcal pharyngitis in preschool children aged 3 months to 5 years. Clinical Pediatrics (June 1999), 38: 357-360.
7. Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA. Group A Streptococcal Pharyngitis in Adults 30 to 65 years of age. Southern Medical Journal (May 1999), 491- 492

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15. BIBLIOGRAFIA

1. Murray, P.R., et al. Manual of Clinical Microbiology, 6th Edition, ASM Press, Washington D.C., 1995, p. 299-307.
2. Webb, K.H. Does Culture Confirmation of High-sensitivity Rapid Streptococcal Tests Make Sense? A Medical Decision Analysis. Pediatrics (Feb 1998), 101:2, 2.
3. Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. Diagnosis and Management of Group A Streptococcal Pharyngitis. Clinical Infectious Diseases (1997), 25: 574-83.
4. Needham CA, McPherson KA, Webb KH. Streptococcal Pharyngitis: Impact of a High-sensitivity Antigen Test on Physician Outcome. Journal of Clinical Microbiology (Dec 1998), 36: 3468-3473.
5. Shea, Y.R., Specimen Collection and Transport, Clinical Microbiology Procedures Handbook, Isenberg, H.D., American Society of Microbiology, Washington D.C., 1.1.1- 1.1.30, 1992.
6. Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. Group A beta-hemolytic streptococcal pharyngitis in preschool children aged 3 months to 5 years. Clinical Pediatrics (June 1999), 38: 357-360.
7. Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA. Group A Streptococcal Pharyngitis in Adults 30 to 65 years of age. Southern Medical Journal (May 1999), 491- 492

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1. Murray, P.R., et al. Manual of Clinical Microbiology, 6th Edition, ASM Press, Washington D.C., 1995, p. 299-307.
2. Webb, KH. Does Culture Confirmation of High-sensitivity Rapid Streptococcal Tests Make Sense? A Medical Decision Analysis. Pediatrics (Feb 1998), 101:2, 2.
3. Bisno AL, Gerber MA, Gwaltney JM, Kaplan EL, Schwartz RH. Diagnosis and Management of Group A Streptococcal Pharyngitis. Clinical Infectious Diseases (1997), 25: 574-83.
4. Needham CA, McPherson KA, Webb KH. Streptococcal Pharyngitis: Impact of a High-sensitivity Antigen Test on Physician Outcome. Journal of Clinical Microbiology (Dec 1998), 36: 3468-3473.
5. Shea, Y.R., Specimen Collection and Transport, Clinical Microbiology Procedures Handbook, Isenberg, H.D., American Society of Microbiology, Washington D.C., 1.1.1- 1.1.30, 1992.
6. Nussinovitch, M, Finkelstein Y, Amir J, Varsano, I. Group A beta-hemolytic streptococcal pharyngitis in preschool children aged 3 months to 5 years. Clinical Pediatrics (June 1999), 38: 357-360.
7. Woods WA, Carter CT, Stack M, Connors Jr AF, Schlager TA. Group A Streptococcal Pharyngitis in Adults 30 to 65 years of age. Southern Medical Journal (May 1999), 491- 492

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