AMNIOQUICK

REF 1090003 / 1090003_C

FR (FRANÇAIS / FRENCH) 1

EN (ENGLISH) 5

DE (DEUTSCH / GERMAN) 9

ES (ESPAÑOL / SPANISH) 13

IT (ITALIANO / ITALIAN) 17

PT (PORTUGUES / PORTUGUESE) 21

CONTENTS

MANUFACTURER INFORMATION

BIOSYNEX S.A.
22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN – France
Standard : **+33 3 88 78 78 87**
www.biosynex.com

Contact France : **+33 3 88 77 57 00**
service.clients@biosynex.com
Contact other countries : **+33 3 88 77 57 52**
sales@biosynex.com
After-sales service: **+33 3 88 77 57 25**
tech.support@biosynex.com

CH REP

BIOSYNEX SWISS S.A.
Route de Rossemaison 100
2800 DELEMONT - Switzerland

NON-STANDARD SYMBOLS

	FR	EN	DE	ES	IT	PT
DIL	Diluant	Diluent	Diluent	Diluyente	Diluente	Diluente
CH REP	Mandataire suisse	Authorized representative in Switzerland	Bevollmächtigter Vertreter in der Schweiz	Representante autorizado en Suiza	Rappresentante autorizzato in Svizzera	Representante autorizado na Suíça

CONTENTS

FR

AMNIOQUICK

TEST DE DETECTION RAPIDE DE L'IGFBP-1 (INSULIN-LIKE GROWTH FACTOR-BINDING PROTEIN 1)

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EN

AMNIOQUICK

RAPID TEST FOR DETECTION OF IGFBP-1 (INSULIN-LIKE GROWTH FACTOR-BINDING PROTEIN 1)

In vitro diagnostic medical device for professional use.

REF 1090003 / 1090003_C

1. INTENDED PURPOSE

AMNIOQUICK is a simple and rapid immunochromatographic test that detects IGFBP-1 (Insulin-like Growth Factor-Binding Protein 1) in vitro based on a vaginal swab. The AMNIOQUICK test is designed to detect foetal membrane rupture in pregnant women using a vaginal swab sample. The test gives a qualitative result.
The AMNIOQUICK test can be used by health professionals in the event of suspected premature rupture of membranes (PROM). Every test device is intended for professional in vitro diagnostic use only.

2. CLINICAL SUMMARY

Premature Rupture of Membranes (PROM) is relatively common and occurs in 5 to 10% of pregnancies. It can lead to premature births and foetal infections. Amniotic fluid leakage cannot always be detected through a physical examination and confirmatory blood tests are often useful. The blood tests involve detecting, in the vagina, either alkaline secretions (a simple and low-cost method, but insufficiently sensitive and specific) or the presence of proteins that can be found in high concentrations in amniotic fluid (diamine-oxidase, alpha-fetoprotein, foetal fibronectin, IGFBP-1).

3. TEST PRINCIPLE

A pair of anti-IGFBP-1 monoclonal antibodies is used to detect IGFBP-1. One is immobilized on the nitrocellulose membrane in the test area and acts as a capture antibody. The other antibody is first coupled with colloidal gold nanoparticles and acts as the revelation reagent. When a sample positive for IGFBP-1 migrates, antigen-antibody complexes form in the sample pad and migrate along the membrane to be captured by the antibodies in the test area, thereby creating a visible purple line in the test area on account of the gold nanoparticles. A purple line in the control area means that the test results are valid and that the procedure has been followed correctly.

4. REQUIRED MATERIAL

Material provided:

REF 1090003	REF 1090003_C
- 10 test strips, individually packaged in an aluminium pouch with a desiccant - 10 sterile Dacron vaginal swabs - 10 diluent, in single-dose vials - 10 patient forms - 1 workstation - 1 instructions for use.	- 10 test strips, individually packaged in an aluminium pouch with a desiccant - 10 sterile Dacron vaginal swabs - 10 diluent, in single-dose vials - 10 patient forms - 1 tube containing a lyophilized bead of positive control - 1 workstation - 1 instructions for use.

Material required but not provided:
Timer

5. PRECAUTIONS

- For in vitro diagnostic use only.
- For best results, carefully follow the procedure and the storage instructions.
- Do not open the aluminium packaging until it is at room temperature to prevent condensation. Damp and high temperatures may affect results.
- Do not use the kit after the expiry date.
- Do not eat, drink or smoke while handling the samples and the test.
- Wear a lab coat, disposable gloves, protective glasses when handling potentially infectious materials and when taking the test.
- The samples and the equipment used when performing the test must be considered as potentially infectious and treated as such. Dispose of the test components and the swab according to the procedure applicable to potentially infectious waste.
- Avoid splatters and formation of aerosols. Clean all spilled liquid with an appropriate disinfectant.
- The tests and tubes provided are intended for single use. Do not re-use the strips, swabs or diluent tubes.
- Do not combine or mix up the reagents of kits or different batches.
- Do not use the test or the control if the aluminium packaging is open or damaged.
- Carefully read the AMNIOQUICK instructions for use before performing the test.
- Do not use a yellow-colored diluent.

6. STORAGE AND STABILITY

Store the kit between +2°C and +30°C. Do not freeze. The kit and its components can be used until the expiry date indicated on the kit label. Do not use after the expiry date. The test strip and positive control must remain in their closed aluminium pouches until use. Just before use, open the pouch and take the test strip or positive control out. After opening, the test strip must be used immediately and the positive control within one hour if kept in its closed Eppendorf tube. Once prepared, the positive control can be kept up to 24

hours at room temperature (+19°C to +24°C).

7. SPECIMENS/CONTROLS COLLECTION AND HANDLING

Collection without speculum:

Use the sterile swab provided to collect secretions on the surface of the vagina. Remove the swab from its packaging. Insert the swab into the vagina (5 cm deep) and soak for 1 minute.

Collection with speculum

Alternatively, a speculum can be used. Collect vaginal secretions after 15 seconds of contact with the vaginal walls in the posterior fornix.

Sample storage

The swab must be diluted in the extraction buffer tube immediately after sampling. The tube can then be stored for 4 hours at most at room temperature or 12 hours at -20°C before performing the test, given that proteases found in vaginal secretions may degrade IGFBP-1.

8. TEST PROCEDURE

Test procedure for raw samples:

1. Bring the complete kit and the samples to room temperature (19-24°C) before performing the test.
2. Open the single-dose vial and place it on a horizontal and flat surface.
3. Immerse the swab into the vial and rotate for a dozen or so seconds. Press the swab against the sides of the tube to extract the most liquid possible from the swab. Remove the swab.
4. Open the aluminium pouch, remove the strip from its packaging and dip it in the vial with arrows pointing downwards. Gently hit the strip on the bottom of the vial to enhance the migration.
5. Start the stopwatch and read the result at 10 minutes. A positive result can appear in the first minutes. No test line at 10 minutes means that the result is negative. Do not interpret after 15 minutes.
6. Dispose of the components of the test and the swab according to the procedure applicable to potentially infectious waste.


Test procedure for positive control (REF 1090003_C):

1. Bring one plastic packet containing components for one test and the pouch containing the positive control bead tube to room temperature (+19°C to +24°C) before performing the test.
2. Open the diluent tube and place it on a horizontal and flat surface.
3. Open the pouch containing the positive control bead tube.
4. Open the tube and drop the positive control bead into the diluent tube.
5. Close the tube and mix by inverting the tube 7 times.
6. Open the aluminium pouch, remove the test strip from its packaging and dip the strip immediately in the tube with the arrow pointing downwards. Gently hit the strip on the bottom of the tube to enhance the migration.
7. Start the timer and read the results at 10 minutes. A positive result can appear in the first minutes. Do not interpret after 15 minutes.
8. Dispose of test strip, diluent tube and positive control tube according to the procedure applicable to potentially infectious waste.


9. RESULTS INTERPRETATION

Results are interpreted as negative, positive or invalid depending on the coloured lines that appear in the control area and test area:


- **POSITIVE:** Two distinct purple lines: a control line appears in the control area and a purple line (even of weak intensity) appears in the test area.

- **NEGATIVE:** A single purple line appears in the control area. No line appears in the test area.

- **INVALID:** No purple line in the control area, regardless of the configuration in the test area. The results of a test without a control line must be discarded. Review the procedure and repeat it with a new test device. If the problem persists, contact the local distributor.

NOTE: Any shade of color in the test and control lines areas should be considered positive. Do not interpret a white "ghost" test line and consider the result as invalid.

10. QUALITY CONTROL

Procedural control:

Internal procedural control is included in the test. A coloured line appearing in the control zone (C) ensures that sufficient specimen volume has been loaded and that the correct procedure has been followed by the operator.

External control:

Good Laboratory Practices recommend the use of control materials to ensure proper kit performance. A positive control is provided within the kit of reference 1090003_C. Follow the test procedure for positive control as described in part 8. According to Good Laboratory Practice, it is recommended to carry out an inspection on each new batch or delivery to confirm the procedure and check performance. Each laboratory must define its own quality control system.

11. PERFORMANCES

Detection limit:

The detection limit for the AMNIOQUICK test is based on a native IGFBP-1 control prepared in the dilution buffer is **2.5 ng/mL**

Clinical studies:

Test performance was determined based on clinical studies. Studies including more than 100 patients are shown in the below table.

	Study n°1	*Study n°2*	*Study n°3*
	Li Yanfang *et al.*	**Newcomb *et al.***	**Wang *et al.***
Place	Beijing, China	Rochester, United States	Chengdu, China
Year	2010	2011	2012
No. of patients	200	272	100
Sensitivity	98.5%	92.2%	100%*
Specificity	100%	93.1%	100%
PPV	100%	92.2%	100%
PNV	97.1%	93.1%	100%*

*Before 37 weeks of amenorrhea

Hook effect:

The AMNIOQUICK test does not have a hook effect for concentrations up to 100 µg/mL of native IGFBP-1 diluted in the extraction tube. However, it is expected that a greater concentration of IGFBP-1 may induce a hook effect.

Reproducibility:

Inter-day:

An inter-day variation study was performed using seven concentrations of a native IGFBP-1 control diluted in the test's diluent on one batch over five days. Results of this study showed that AMNIOQUICK reproducibility between days with tested negative and positive samples is 100%.

Inter-batch:

An inter-batch variation study was performed using seven concentrations of a native IGFBP-1 control diluted in the test's diluent on three batches. Results of this study showed that AMNIOQUICK reproducibility between batches with tested negative and positive samples is 100%.

Inter-operator and inter-place:

An inter-operator variation study was performed using seven concentrations of a native IGFBP-1 control diluted in the test's diluent on one batch by two different operators at two different places. Results of this study showed that AMNIOQUICK reproducibility between operators and places with tested negative and positive samples is 100%.

Interference:

The following substances were added to the diluent (0.5mL), then tested:

- Blood (20 and 50 µL)
- Serum (50 µL)
- Urine of pregnant women (50 µL)
- Sperm (50 µL)
- Lubricants (50 µL)
- C vitamin (1.25 g)

This study showed that the presence of a high quantity of urine (100 µL/mL) does not impact the AMNIOQUICK test results.
However, the presence of 100 µL/mL of serum can lead to false-positive results. Moreover, the presence of 40 µL/mL of blood creates a background noise that makes results difficult to read. For those two reasons, results obtained with bloody samples must be considered with caution. The lubricants (K-Y®, KLY® and CerviLube™) tested at concentration of 100 µL/mL do not interfere with AMNIOQUICK test results. Sperm does not lead to any interference up to 100 µL/mL.
The presence of up to 2.5 g/mL of C vitamin does not impact the test results.

12. LIMITATIONS

1. As for any diagnostic test, the test result must be correlated with clinical results.
2. If there is a lot of blood (head of the swab is noticeably red) in the vaginal secretion sample, a positive result must be interpreted with caution.
3. False negative results may occur if the test is performed 12 hours after stopping amniotic fluid leakage due to membrane rupture.
4. Collecting a sample incorrectly may lead to false negative results.
5. In case of frank membrane rupture, the use of AMNIOQUICK is not recommended.
6. AMNIOQUICK is a qualitative test. No quantitative interpretation can be made based on the test result.

13. LITERATURE

1. Friedman ML, McElin TW. Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol 1969:104:544-550
2. Steinman G, Kleiner GJ, Greston WM. Spontaneous rupture of membranes. NY St. J Med 1979:1849-1851
3. Koninckx PR, Trappeniers H, van Assche FA. Prolactin concentration in vaginalfluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol 1981;88:607–610
4. Rochelson BL, Richardson DA, Macri JN. Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol 1983;62:414–418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol 1986:93:909-915
6. Suikkari AM, Rutanen EM, Seppällä M. Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod 1987:2:297-300
7. Bell SC, Keyte JW; N-terminalo amnio acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. The Endoc Soc 1988:1202-1204
8. Rutanen EM, Pekonen F, Kärkkainen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vagonal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta 1993:214:73-81
9. Medina TM, Hill DA. Preterm premature rupture of membranes: diagnosis and management. Am Famil Physic 2006:73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol 1994;171:146–150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. Détection rapide de l'Insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem 1999 :36 :388-390
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. Insulin-Like Growth Factor Binding Protein-I; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci 2005:35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin 2005 :12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. The value of insulin-like growth factor binding protein-1 in the cervical-vagonal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol 2005:121:159-163
15. CNGOF (Collège national des gynécologues et obstétriciens français). Recommandations pour la partie clinique. Rupture prématurée des membranes avant terme. 2018. Gynécologie Obstétrique Fertilité & Sénologie. https://doi.org/10.1016/j.gofs.2018.10.016

14. MODIFICATION HISTORY

IFU version	Amended paragraph	Details of the changes
V3	Content	Addition of the table of contents and reference 1090003_C.
	4 : Required material	Addition of lyophilized positive control in ref 1090003_C.
	5 : Precautions	Addition of a precaution.
	6 : Storage and stability	Full review of this section.
	7 : Specimens/controls collection and storage	Update in sample storage data.
	8 : Test procedure	Addition of procedure for positive control.
	9 : Results interpretation	Addition of a note.
	10 : Quality control	Addition of the positive control within the kit of reference 1090003_C. Removal of ref 6090001.
	11 : Performances	Update of detection limit. Update. Correction in clinical study 3 data. Update on hook effect and interfering substances. Update on reproducibility studies.

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