CHLAMYTOP

REF 1030003

CE 0123

FR (FRANÇAIS / FRENCH) 1

EN (ENGLISH) 6

DE (DEUTSCH / GERMAN) 11

ES (ESPAÑOL / SPANISH) 17

IT (ITALIANO / ITALIAN) 23

PT (PORTUGUES / PORTUGUESE) 29

CS (ČESKY / CZECH) 35

HU (MAGYAR / HUNGARIAN) 40

PL (POLSKI / POLISH) 46

SK (SLOVENSKA / SLOVAK) 52

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EN

CHLAMYTOP

RAPID TEST FOR THE QUALITATIVE DETECTION OF CHLAMYDIA ANTIGENS.

In vitro diagnostic medical device for professional use.

REF 1030003

1. INTENDED PURPOSE

CHLAMYTOP is an in-vitro diagnostic lateral flow immunoassay intended for the qualitative detection of Chlamydia antigens from female endocervical swabs and male urethral swabs. This rapid test is an aid in the diagnosis of chlamydia infection.

2. CLINICAL SUMMARY

The genus Chlamydia is made up of three species: Chlamydia trachomatis, Chlamydia pneumoniae (1), a primarily human pathogen, and Chlamydia psittaci, a primarily animal pathogen. Chlamydia trachomatis consists of 15 known serovars. These are associated with urogenital infections and lymphogranuloma venereum (LGV). Chlamydia trachomatis infections are the most common bacterial sexually transmitted diseases. In the USA there are about 4 million new cases every year, mainly cervicitis and non-gonococcal urethritis (8). This organism also causes conjunctivitis and pneumonia in children (2, 4–7). Chlamydia trachomatis infections have both a high prevalence and asymptomatic carriage rate, with frequent serious complications in both women and neonates. Complications of chlamydia infections in women include cervicitis, urethritis, endometritis, pelvic inflammatory disease (PID) and increased incidence of ectopic pregnancy and infertility. Vertical transmission of the disease from the mother to the neonate during parturition can result in inclusion conjunctivitis and pneumonia (8). In men at least 40% of all cases of non-gonococcal urethritis are associated with Chlamydia infections and epididymitis (6). Approximately 70% of women with endocervical infections and up to 50% of men with urethral infections are symptomatic (8). Chlamydia psittaci infections are associated with respiratory diseases in individuals exposed to infected birds and are not transmitted from human to human. Chlamydia pneumoniae, which was first isolated in 1983, is associated with respiratory infections and pneumonia (2). Traditionally chlamydia infections have been diagnosed by the detection of chlamydiae in tissue culture cells or polymerase chain reaction (PCR). Together with PCR, the culture method is the most sensitive and specific laboratory method, but it is labour-intensive and expensive, it takes a long time (2–3 days) and is not routinely available in most institutions (2, 3, 7). Direct tests such as immunofluorescence assay (IFA) require special equipment and a skilled operator to read the results.

3. TEST PRINCIPLE

CHLAMYTOP rapid test has been designed to detect *Chlamydia* through visual interpretation of color development in the internal strip. The membrane was immobilized with monoclonal antibodies to *Chlamydia trachomatis* lipopolysaccharide (LPS) on the test region (T). During the test, the specimen is allowed to react with latex-colored monoclonal anti-Chlamydia antibodies, which were precoated on the conjugate pad of the test. The mixture then moves on the membrane by capillary action, and interact with reagents on the membrane. If there were enough Chlamydia antigens in specimens, a colored band will form at the test region (T) of the membrane. Appearance of a colored band at the control region (C), which is coated with proteins, serves as a procedural control: this indicates that proper volume of specimen has been added and membrane wicking has occurred. Buffer A serves to crack the sample so that the antigen can be fully bound. Buffer B serves to neutralize the pH so that the sample is in the right environment.

4. REQUIRED MATERIAL

Material provided:

- Test cassettes individually packaged
- Extraction solution A (0.2M NaOH) - white cap

Warning
Causes skin irritation.
Causes serious eye irritation.

- Extraction solution B (0.2M hydrochloric acid) - yellow cap
- Extraction tubes and dropper caps
- Female swab: sterile swabs with Dacron tips

	EU REP	
Heinz Herenz Medizinalbedarf GmbH Rudorffweg 10, 21031 Hamburg - GERMANY	/	CE 0482
Jiangsu Rongye Technology Co., Ltd. Touqiao Town, Yangzhou City 225109 Jiangsu, P.R. China	Riomavix S.L. Calle de Almansa 55, 1D, Madrid 28039, Spain	CE 0197

- Reagent holder

Material required but not provided:

- Timer
- Male swab: sterile standard wire-shafted fiber-tipped swab

5. PRECAUTIONS

- For single use only. Do not reuse test cassettes, swabs, extraction tubes and dropper caps.
- Do not use any kit components after the expiry date stated on the kit label.
- Do not use test or swab if pouch is damaged.
- Do not use if there is evidence of microbial contamination or precipitation. Biological contamination can lead to false results.
- Do not mix kit components from different lots. Do not mix extraction solution bottle caps. Do not separate the test kit components. Do not disassemble the cassette. Put back extraction solutions in the kit after each use.
- Samples and used materials must be considered as potentially infectious and handled as such. Dispose of used components according to the procedure and local regulations applicable to potentially infectious waste. Clean surface and any spilled liquid with an appropriate disinfectant.
- Do not eat, drink or smoke in the area where specimens and kit reagents are handled. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when collecting and testing specimens.
- Humidity and temperature can adversely affect results.
- Extraction solution A: Wash hands thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment.
- For endocervical sampling, use only Dacron or Rayon tipped sterile swabs with plastic shafts. It is recommended to use the swab supplied within the kit. Other swabs have not been validated. Swabs with cotton tips or wooden shafts are not recommended. Do not use cytology brushes with pregnant patients.
- For urethral sampling, standard wire-shafted fiber-tipped swabs should be used. The swab for urethral sampling is not provided within the kit.
- Safety data sheet is available upon request.

6. KIT STORAGE

The test kit should be stored at +2°C to +30°C until the expiry date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE.

7. SAMPLE/CONTROLS COLLECTION AND STORAGE

The quality of specimen obtained is of extreme importance. Detection of chlamydia requires a rigorous and thorough collection technique which provides cellular material rather than just body fluids.

For cervical specimens:

1. Use the swab provided within the kit.
2. Before specimen collection, remove excess mucus from the endocervical area with a separate swab or cotton ball and discard. The swab should be inserted into the cervical canal, past the squamocolumnar junction, until most of the tip is no longer visible. This will permit acquisition of columnar or cuboidal epithelial cells which are the main reservoir of the Chlamydia organisms. Firmly rotate the swab for 15 to 20 seconds without contaminating it with exocervical or vaginal cells.
3. Put the swab into the extraction tube if the test is to be carried out immediately after the specimen is collected.

For urethral specimens:

1. Standard wire-shafted fiber-tipped swabs should be used for urethral specimen collection. Instruct the patients not to urinate for at least two hours prior to specimen collection.
2. Insert the swab approximately 2 to 4 cm into the urethra, rotate for 3 to 5 seconds and remove it. Put the swab into the extraction tube if the test is to be carried out immediately after the specimen is collected.

Storage of cervical and urethral samples:

1. Do not place the swab in any transport device containing medium since transport medium interferes with the test and viability of the organisms is not required for this test. If immediate testing is not possible, the patient samples should be placed in a dry transport tube for storage or transport.
2. The swabs can be stored for 24 hours at room temperature (+15 to +30°C) or 1 week at +4°C or up to 6 months at -20°C. All samples should be at room temperature (+15 to +30°C) before the test is performed.

8. TEST PROCEDURE

Bring kit components and samples to room temperature (+15 to +30°C) before use.


Sample extraction steps:

1. Place a clean and identified extraction tube into the holder. Add 8 drops (approx. 400µl) of extraction solution A to the extraction tube.
2. Immerse the patient' swab in the extraction tube and extract for 2 minutes at room temperature. During extraction, use a circular motion to roll the swab against the side of the extraction tube so that the liquid is expressed from the swab and can be reabsorbed.
3. At the end of the extraction time, add 8 drops (approx. 400µl) of the extraction solution B to the extraction tube and extract for another 1 minute in the same way. Then press the swab firmly against the tube in order to expel as much liquid as possible from the swab. Discard the swab in accordance with guidelines for handling infectious materials.
4. Put the dropper cap on the extraction tube.
5. The specimens extracted can be kept at room temperature for 60 minutes without affecting the test results.

Test handling steps:


1. Remove the test cassette from its sealed pouch and place it on a clean, dry and level surface. Label the test cassette with patient identification. The test must be performed within one hour.
2. Add 3 drops (approx. 100µl) of the extracted sample from the extraction tube to the sample well (S) of the test cassette. Avoid trapping air bubbles in the sample well (S), and do not drop any solution in observation window.
3. Wait for the colored line(s) to appear. The test result should be read at 15 minutes. Do not interpret the result after 20 minutes.

9. INTERPRETATION OF RESULTS


POSITIVE:
Two colored lines appear. One line appears in the control line area (C) and one line in the test line area (T). A positive result indicates that chlamydia antigen has been detected.

NOTE: The intensity of colour in the test area (T) may vary depending on the concentration of chlamydia antigens present in the specimen. Therefore, any shade of colour in the test area (T) should be considered positive.

NEGATIVE:
One colored line appears in the control line area (C). No colored line appears in the test line area (T). A negative result indicates that no chlamydia antigen is present in the specimen or that it is below the detection level of the test.

INVALID:
Control line (C) fails to appear. Insufficient specimen volume, expired test components or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the kit immediately and contact your local distributor.

10. QUALITY CONTROL

Internal controls:

A procedural internal control is included in the test. A colored line appearing at the control line (C) is considered a positive procedural internal control, confirming sufficient sample volume and correct procedure.

External controls:

Good laboratory practice includes the use of external controls to ensure proper kit performance. It is recommended to perform a positive and negative external control for every kit, and as deemed necessary by internal laboratory procedures. External positive and negative controls are not provided with the kit. To perform a control test, complete the steps in the Test Procedure section (8 treating the control swab in the same manner as a specimen swab.

11. PERFORMANCES

Clinical sensitivity and specificity:

- *Comparison with another Rapid Diagnostic Test*

A clinical study compared results of the CHLAMYTOP test to results obtained with another Rapid Diagnostic Test. Discrepant results were evaluated with Real-Time PCR. 558 cervical specimens and 528 urethral specimens were tested.

Cervical specimens:

		Other Rapid Test (+ PCR)		
		Positive	Negative	Total
CHLAMYTOP	Positive	80	4	84
	Negative	5	469	474
	Total	85	473	558

Relative Sensitivity = 94.12% (95% CI: 86.80%-98.06%)*
Relative Specificity = 99.15% (95% CI: 97.85%-99.77%)*
Overall Agreement = 98.39% (95% CI: 96.96%-99.26%)*
**CI: Confidence Interval*

Urethral specimens:

		Other Rapid Test (+ PCR)		
		Positive	Negative	Total
CHLAMYTOP	Positive	105	0	105
	Negative	6	417	423
	Total	111	417	528

Relative Sensitivity = 94.59% (95% CI: 88.62%-97.99%)*
Relative Specificity ≥ 99.9% (95% CI: 99.12%-100%)*
Overall Agreement = 98.86% (95% CI: 97.54%-99.58%)*
**CI: Confidence Interval*

- *Comparison with a CE-marked PCR test*

A clinical study compared results of the CHLAMYTOP test to results obtained with a CE-marked PCR test. Discrepant results were evaluated with another Real-Time PCR test. 193 cervical specimens and 162 urethral specimens were tested

Cervical specimens:

		PCR reference method		
		Positive	Negative	Total
CHLAMYTOP	Positive	51	3	54
	Negative	13	126	139
	Total	64	129	193

Relative Sensitivity = 79.69% (95% CI: 67.77%-88.72%)*
Relative Specificity = 97.67% (95% CI: 93.35%-99.52%)*
Overall Agreement = 91.71% (95% CI: 86.89%-95.19%)*
*CI: Confidence Interval

Urethral specimens:

		PCR reference method		
		Positive	Negative	Total
CHLAMYTOP	Positive	51	1	52
	Negative	10	100	110
	Total	61	101	163

Relative Sensitivity = 83.61% (95% CI: 71.91%-91.85%)*
Relative Specificity = 99.01% (95% CI: 94.61%-99.97%)*
Overall Agreement = 93.21% (95% CI: 88.18%-96.56%)*
**CI: Confidence Interval*

Analytical sensitivity:
The detection level of CHLAMYTOP is 1.0 x 10^5 IFU/mL, based on spiked samples.

Precision:
Three batches of CHLAMYTOP has been tested and give repeatable results. As well, inter days, lots, places and persons reproducibility has been validated on three batches.

Interferences:
No interference has been observed between CHLAMYTOP test results and following substances:

Metronidazole suppositories	10mg/ml Human mucin
Miconazole Nitrate suppositories	Tindazole vaginal effervescent tablets
Confort Pessaries	Bifonazole cream
Nonoxynol pellicles Dissolvable spermicidal paper	Mycostatin Effervescent vaginal tablets nystatin
Fuyinjie Lotion	Jieeryinxiye Lotion
Lubricating jelly	Disinfector

An interference can be observed between CHLAMYTOP test and human whole blood at a concentration >40µL in case of female swab and >10µL in male swab (see Limitations).

Cross reactions:
The antibody used in CHLAMYTOP is targeted at a genus-specific epitope, which is present in all 15 Chlamydia serovars. In addition, *Chlamydia psittaci* and *Chlamydia pneumoniae* strains have been tested with CHLAMYTOP and gave positive results.
Cross reactivity with other organisms has been studied using suspensions of 10^7 IFU/mL specimen material. The following organisms were not detected with CHLAMYTOP:

Acinetobacter calcoaceticus	*Acinetobacter spp.*	*Streptococcus faecalis*
Proteus vulgaris	*Gardnerella vaginalis*	*Neisseria catarrhalis*
Neisseria lactamica	*Staphylococcus aureus*	*Pseudomonas aeruginosa*
Salmonella typhi	*Streptococcus faecium*	*Neisseria gonorrhoea*
Escherichia coli	*Candida albicans*	*Trichomonas vaginalis*
Ureaplasma urealyticum	*Mycoplasma hominis*	*Neisseria meningitidis*

Hook effect:
No hook effect has been observed up to a concentration of 10^9 IFU/ml.

12. LIMITATIONS

- The test should be used for the qualitative detection of Chlamydia antigens from female endocervical swabs and male urethral swabs. The chlamydia antigens concentration cannot be determined by this test.
- The test does not specifically differentiate between *C.trachomatis, C.pneumonia or C.psittaci*.
- Detection of chlamydia is dependent on the concentration of organisms present in the specimen. This may be affected by specimen collection methods and patient factors such as age, previous history of STDs, presence of symptoms, etc. The minimum detection level of this test may vary according to serovar.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on one single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- A negative result does not exclude chlamydia infection when the number of infectious agents is below the detection limit.
- Excessive blood in sample (>40µL in case of female swab and >10µL in male swab) may cause false positive results. Endocervical samples from female patients should not be collected during menstrual period. If the swab is collected with excessive blood or clots, it needs to be resampled.

13. LITERATURE

1. Grayston, J.T. et al., N. Engl. J. Med. 315:161, 1986.
2. Barnes, R., Clin. Microbiol. Rev. 2:119, 1989.
3. Kellogg, J.A., Arch.Pathol. Lab. Med. 113:453, 1989.
4. Schachter, J., N., Engl. J. Med. 13:105, 1978.
5. Schachter, J., Chlamydiae, p.1045ff, Manual of Clinical Microbiology, 5th ed, ASM, Washington,1991.
6. Schachter, J. & Dawson, C.R., Sex. Transm. Dis. 8:167, 1981.
7. Stamm, W.E., Ann. Intern. Med. 108:710, 1988.
8. CDC. Morbid.Mortal. Weekly Rep. 34(suppl 3-s): 535, 1985

14. MODIFICATION HISTORY

IFU version	Amended paragraph	Details of the changes
V1	All	Multilingual format and transition to Version 1 Update IFU to Biosynex format
	Cover page	Addition of the table of contents
	4: Required materials	Replacement of the EC Rep symbol with the EU Rep symbol for the swab
	14: Modification history	Addition of this paragraph
V2	4: Required materials	Removal of the instructions for use

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