# STREPTATEST

REF 10025ND

FR (FRANÇAIS / FRENCH) ..... 1

EN (ENGLISH) ..... 7

DE (DEUTSCH / GERMAN) ..... 13

IT (ITALIANO / ITALIAN) ..... 19

NL (NEDERLANDS / DUTCH) ..... 25

NO (NORSK / NORWEGIAN) ..... 31

DA (DANSK / DANISH) ..... 36

FI (SUOMALAINEN / FINNISH) ..... 41

SV (SVENSKA / SWEDISH) ..... 46

CONTENTS

## MANUFACTURER INFORMATION

**BIOSYNEX S.A.**
22 boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN – France
Standard : **+33 3 88 78 78 87**
www.biosynex.com

Contact France : **+33 3 88 77 57 00**
service.clients@biosynex.com
Contact other countries : **+33 3 88 77 57 52**
sales@biosynex.com
After-sales service: **+33 3 88 77 57 25**
tech.support@biosynex.com

CH REP **BIOSYNEX SWISS S.A.**
Route de Rossemaison 100
2800 DELEMONT - Switzerland

## NON-STANDARD SYMBOLS

| | FR | EN | DE | IT | NL | NO | DA | FI | SV |
|---|---|---|---|---|---|---|---|---|---|
| DIL | Diluant | Extraction buffer | Verdünnungspuffer | Tampone di diluizione | Extractiebuffer | Ekstraksjon buffer | Buffer til ekstraktion | Uuttopuskuri | Buffert för extrahering |

## CONTENTS

FR

# STREPTATEST

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## 14. BIBLIOGRAPHIE

1. American Academy of Pediatrics. Peter, G., éd. 1994 Red Book: Report of the Committee on Infectious Diseases. 23rd éd. Elk drove Village, IL; American Academy of Pediatrics; 1994: p. 433.
2. Bisno AL. Group A streptococcal infections and acute rheumatic fever.N. Engl. J. Med. 325: 783-793(1991).
3. Bisno AL, Pearce IA, Wall HP, Moody MD, and Stollerman GH. Contrasting epidemiology of acute rheumatic feverand acute glomerulonephritis. N. Eng. J. Med. 283: 561-565 (1970).
4. Edwards, E.A., Phillips, LA., and Suiter, W.C., Diagnosis of Group A Streptococcal Infections Directly from Throat Sécrétions, J. Clin. Micro., 15,481-483(1982).
5. Facklam RR. U.S. Dept. of Health and Human Services, PHS, CDC, Pub. No. CDC 77-13.
6. Gupta, R., Talwar, G.P. and Gupta, S.K., Rapid Antibody Capture Assay for Détection of Group A Streptococci Using Monoclonal Antibody and Colloïdal Gold Monospecific Polyvalent Antibody Conjugale, J. Immunoassay, 13,441-445(1992).
7. Heiter, BJ, Bourbeau PP, Comparison of two rapid streptococcal antigen détection assays with culture for diagnosis of streptococcal pharyngitis. J. clin. Microb. 16 : 748-53 (1995)
8. Kuttner AG and Krumwiede E. Observations on thé effect of Streptococcal upper repiratory infections on rheumatic children. a three-year study. J. Clin. Invest. 20: 273-287 (1941).
9. Lauer BA, Relier LB and Mirrell S. Effect of atmosphère and duration of incubation on primary isolation of group A Streptococci from throat cultures. J. Clin. Microb. 17:338-340 (1983)
10. Levinson, M.L. and Frank, PP., Differentiation of Group A from other Bêta Hemolytic Streptococci with Baci-tracin, J. BacterioL 69, 284-287 (1955).
11. Potter EV, Svartman M, Mohamed 1, Cox R, Poo-King T, and Earle DP. Tropical acute rheumatic fever and associated streptococcal infections compared with concurrent acute glomerulonephritis. J. Pediatr. 92: 325-333 (1978).
12. Ruoff KL, Whiley RA, Beighton D. Streptococcus. In: Manual of Clinical Microbiology, 7th éd., PR Muray, EJ Baron (eds), American Society of Microbiology, Chapter 17, pp. 283-297 (1999).
13. Ross, P.W., Throat Swabs and Swabbling Technique, The Practitioner, 207, 791-796(1971).
14. Wannamaker LW. Différences between streptococcal infections of thé throat and skin. N. Eng. J. Med. 282:78-85(1970).
15. Wannamaker LW. Changes and changing concepts in thé biology of group A Streptococci and thé epidemiology of streptococcal infections. Rev. Infect. Dis., 2: 967-973, (1979).
16. Poyart C., Bingen E., Maisonneuve P., Deschenes M., Charlier-Bret N., Boucher B. Réévaluation des Tests de Diagnostic Rapide des Angines à Streptocoque du Groupe A (Streptococcus pyogenes) (Dec. 2002). http://afssaps.santé.fr/htm/5/reactif/tdr.pdf.

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EN

# STREPTATEST

RAPID TEST FOR THE QUALITATIVE DETECTION OF GROUP A STREPTOCOCCI FROM OROPHARYNGEAL SWABS

*In vitro diagnostic medical device for professional use.*

REF 10025ND

## 1. INTENDED USE

STREPTATEST is a rapid lateral flow immunoassay for the qualitative, presumptive detection of Group A Streptococcus antigen from oropharyngeal specimens. This kit is intended for use as an aid in the diagnosis of Group A streptococcal infections.

## 2. INTRODUCTION

Group A beta-haemolytic streptococcus is one of the main germs responsible for upper respiratory infections, including throat infections, pharyngitis and scarlatina. It is important to differentiate Group A streptococcal infections from other infections (e.g. viral) in order to initiate an appropriate treatment. Early diagnosis and treatment of Group A streptococcal infections have made it possible to reduce the severity of symptoms and the number of complications such as glomerulonephritis and acute articular rheumatism. Traditional methods require 24 to 48 hours to reveal and identify the germ. Immunochromatography techniques make it now possible to detect the specific Group A streptococcus antigen directly from a swab. This enables the practitioner to make the diagnosis immediately and prescribe the right treatment.

## 3. TEST PRINCIPLE

STREPTATEST is a lateral flow immunoassay on a membrane using a sandwich capture method. An anti-streptococcal antibody for the A strain is fixed on the test region of the membrane. A second anti-streptococcal antibody for the A strain is conjugated with mauve latex particles and placed just above the membrane immersion zone.
First, the specific antigen for Group A streptococci is extracted from the swab using extraction reagents. The lower part of the strip is then immersed in the extraction solution. The specific Group A streptococcus antigen will bind to the antibody marked with latex particles. The mixture will migrate by chromatography along the membrane and the complex will be fixed in the test zone.
The presence of a mauve line in the test zone indicates a positive result whilst the absence of a mauve line indicates a negative result. The appearance of a mauve band in the control zone reveals that the test has functioned correctly. The absence of this band indicates that the test has not functioned correctly, and the result is not valid.

## 4. MATERIALS

**Materials provided**

- Aluminium pouches containing one strip and one desiccant
- Tongue depressors
- Sterile swabs (CE 0197)
- Extraction tubes with dropper caps
- Extraction reagent A (Sodium Nitrite 2 M), 10 mL
- Extraction reagent B (Acetic Acid 0.2 M), 10 mL
- Positive control (inactivated Strep A), 1 mL
- Negative control (inactivated Strep C), 1 mL

**Reagent holder** **Materials required but not provided**

- Timer

## 5. STORAGE AND STABILITY

All the elements contained in the box can be stored at room temperature or in the refrigerator (+2°C to +30°C). DO NOT FREEZE. Do not use beyond the expiry date indicated on the box.

## 6. PRECAUTIONS

- Carefully follow instructions for use.
- For single use only.
- Do not use test or swab if pouch is damaged.
- Samples and used materials must be considered as potentially infectious and treated as such. Observe established precautions against microbiological hazards during the procedure and follow standard procedures for sample disposal. Dispose of used components according to the procedure and local regulations applicable to potentially infectious waste. Clean surface and any

spilled liquid with an appropriate disinfectant.

- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Humidity and high temperature can adversely affect results.
- Do not use more than the required amount of liquid.
- Do not spill the specimens into the reaction area. Do not touch the reaction area of the device to avoid contamination.
- Do not mix reagents from different lots. Do not mix reagent bottle caps.
- Use only the swab provided in the test kit.
- Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.
- Safety datasheet available on demand.
- Labeling and pictograms according to Regulation (EC) No 1907/2006 (REACH):

| Pictogram | Components | Hazard and precautionary statements (CLP): |
|---|---|---|
|  | Reagent A | Warning<br>H302: Harmful if swallowed.<br>P264: Wash hands, forearms and face thoroughly after handling.<br>P501: Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation. |
|  | Negative control<br>Positive control | Warning.<br>H317 - May cause an allergic skin reaction.<br>H412 - Harmful to aquatic life with long lasting effects.<br>P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.<br>P273 - Avoid release to the environment.<br>P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.<br>P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.<br>P362+P364 - Take off contaminated clothing and wash it before reuse.<br>P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation. |

## 7. SPECIMEN COLLECTION AND STORAGE

1. Depress the tongue using a tongue depressor to avoid contaminating the swab with saliva. The tonsils, the pharynx and all the inflammatory areas, whether ulcerative or exudative, must be swabbed.
2. Performing the test as soon as possible is recommended. If not possible, the swabs may be stored for 4 hours at room temperature (15°C to 30°C) in a dry, sterile and hermetically sealed receptacle or for 24 hours in the refrigerator. If a culture is also requested, use a second swab.

## 8. TEST PROCEDURE

Bring tests, extraction solutions, swab specimens, and/or external controls to room temperature (15-30°C) before testing.

1. Remove the strip from its aluminium pouch just before performing the test.
2. Place **4 drops** of the pink-colored Extraction reagent A into the extraction tube, then add **4 drops** of the colorless Extraction reagent B. Mix the two solutions by gentle stirring. The colour of the final mixture must turn from pink to colourless.
3. Introduce the swab into the tube. Stir the swab by rotating it about ten times in the extraction mixture. **Leave the swab** to rest for **one minute** in the extraction tube.
4. Squeeze the swab strongly against the tube walls in order to expel as much liquid as possible. Discard the swab.
5. Immerse the strip in the extraction tube with the arrows pointing towards the extraction solution. Leave the strip in the extraction tube.
6. **Read the result after 5 minutes.** Depending on the germ concentration, a positive result may appear within the first minutes. Nevertheless, 5 minutes are required for a negative result to be confirmed. Do not read beyond 10 minutes.

*Collect the sample using the swab provided*

*4 drops of Reagent A*

*4 drops of Reagent B*

*Rotate the swab 10 times into the tube*

1 min.


*Squeeze the swab before discarding it*

*Immerse the test strip into the tube*

5 min.

## 9. RESULTS INTERPRETATION

Two red-coloured bands can be seen both in the control zone and in the test zone.

*NOTE: The intensity of colour in the test area (T) may vary depending on the concentration of Strep A present in the specimen. Therefore, any shade of colour in the test area (T) should be considered positive.*

Only one red-colored band appears in the control zone.

The absence of any coloured band in the control zone indicates that something went wrong with the test. Review the procedure and repeat the test with new components. If the problem persists, discontinue using the test kit and contact your local distributor.

## 10. QUALITY CONTROL

***Internal controls:***

A procedural internal control is included in the test. The control band is considered as an internal control of the good functioning of the strip and the reagents.

***External controls:***

Testing the positive and negative controls is recommended for each new batch delivery and each new user. Besides, the positive and negative controls should be tested whenever the integrity of the test (storage, procedure) is questioned. External positive and negative controls are supplied with the kit.

Procedure for External Quality Control Testing:

1. Add **4 drops** of Extraction reagent A and **4 drops** of Extraction reagent B into an extraction tube. Mix vigorously.
2. Add **1 drop** of positive or negative control solution into the tube.
3. Introduce a clean swab and follow the test procedure from step 3 of the test procedure part.
4. Positive control is expected to give a positive result. Negative control is expected to give a negative result. Please, contact your local distributor if unexpected control results are obtained.

## 11. LIMITATIONS

1. The quality of the test depends on the quality of the sample collected. A false negative result may be the consequence of poor sample collection or inappropriate swab storage. A negative result may also be obtained in patients at early stage of infection due to insufficient antigen concentration. Therefore, a new sample must be collected and tested by traditional culture methods if a Group A streptococcus infection is suspected despite a negative result.
2. In rare cases, samples heavily colonized with Staphylococcus aureus may give false positive results. If the clinical signs are not in agreement with the test result, performing a traditional bacteriological examination is recommended.
3. Respiratory infections, including pharyngitis, may be due to streptococci belonging to serogroups other than A as well as from other pathogens.
4. STREPTATEST does not allow a quantitative assessment of the concentration of Group A streptococci.
5. As for any *in vitro* diagnosis test, the clinical diagnosis should not be based on the sole test result but should be made by the clinician after all clinical and laboratory evaluations have been completed.

## 12. EXPECTED VALUES

10 to 25% of adult throat infections and 25 to 40% of children throat infections are due to Group A beta-haemolytic streptococci. The infection is particularly prevalent in winter and at the beginning of spring in areas with high urban concentrations.

## 13. PERFORMANCE CHARACTERISTICS

***Diagnostic Sensitivity and Specificity***

A comparative study was performed using 525 patients showing signs of throat infection of pharyngitis. Each sample was seeded on blood agar. Beta-haemolytic colonies obtained from agar culture were then confirmed as being Group A streptococci by grouping techniques. Each swab was also tested with STREPTATEST. 26 patients were excluded from the study because the data collected was incomplete. Results were summarised as follow.

| | | **STREPTATEST** | | Total |
|---|---|---|---|---|
| | | Positive | Negative | |
| **CULTURE** | Positive | 120 | 4 | 124 |
| | Negative | 20 | 355 | 375 |
| | Total | 140 | 359 | 499 |

**Sensitivity 96.8 %[95%CI* : 91 %-99 %]**
**Specificity 94.7 %[95%CI*: 92 %-97 %]**
**PPV 85.7% [95%CI*: 79 %-91 %]**
**NPV 98.9% [95%CI*: 97 %-100 %]**
**Correlation 95.2 %[95%CI*;93 %-97 %]**
***Confidence interval**

Considering the different colonisation levels, the results are distributed as follows:

| Culture classification | | Culture results | |
|---|---|---|---|
| | | STREPTATEST/Culture | Percentage |
| Sensitivity | Rare (< 10 colonies) | 10/11 | 91% |
| | 1+ | 9/9 | 100% |
| | 2+ | 17/19 | 89% |
| | 3+ | 36/37 | 97% |
| | 4+ | 48/48 | 100% |
| | Total sensitivity | 120/124 | 96.8% |
| Specificity | | 355/375 | 94.7% |
| Global correlation | | 475/499 | 95.2% |

***Detection limit***

In an initial study, the detection limit of STREPTATEST was studied on 8 different strains of Group A streptococci. 5 strains presented a detection limit close to $1x10^4$ bacteria/test. The 3 other strains presented a detection limit close to $1x10^5$ bacteria/test (see table below).

| **ATCC STRAIN NUMBER** | **DETECTION LIMIT** |
|---|---|
| Do not reutilize12202 | $2.5 \times 10^5$ bacteria/test |
| 12203 | $3.13 \times 10^4$ bacteria/test |
| 12204 | $3.13 \times 10^4$ bacteria/test |
| 12365 | $6.25 \times 10^4$ bacteria/test |
| 14289 | $3.13 \times 10^4$ bacteria/test |

| | |
|---|---|
| 19615 | 3.13 x $10^4$ bacteria/test |
| 49399 | 5.0 x $10^5$ bacteria/test |
| 51399 | 1.25 x $10^5$ bacteria/test |

The results obtained in a second external study carried out on 4 strains were summarized in the following table:

| STRAIN | STRAIN CHARACTERISTICS | DETECTION LIMIT (UFC/mL) | DETECTION LIMIT (bacteria/test) |
|---|---|---|---|
| ATCC 19615 | Reference strain | $10^5$ | 5x$10^3$ |
| ATCC 12344 (STM1) | Wild type strain | $10^6$ | 5x$10^4$ |
| NCTC 8328 | Scarlatina strain | $10^5$ | 5x$10^3$ |
| ATCC 12358 (STM 19) | Rheumatic strain | $10^5$ | 5x$10^3$ |

Conclusion: The STREPTATEST detection limit ranges from 5.0 x $10^3$ bacteria/test to 5.0 x $10^5$ bacteria/test depending on the studies and the strains analysed.

***Accuracy***

**Intra-test**: The intra-test accuracy was determined on 15 replicates of 3 samples (one negative, one weak positive and one strong positive). The 3 samples were correctly identified in more than 99% of cases.

**Inter-test**: The inter-test accuracy was determined from 15 independent tests on 3 samples (negative, weak positive and strong positive). 3 batches of STREPTATEST were used. The 3 samples were correctly identified in more than 99% of cases.

***Cross Reactivity***

The following organisms were tested at 1.0 x $10^7$ organisms per test and were all found to be negative when tested with STREPTATEST.

| | | |
|---|---|---|
| *Group B Streptococcus* | *Streptococcus sanguis* | *Neisseria sicca* |
| *Group C Streptococcus* | *Serratia marcescens* | *Neisseria subflava* |
| *Group F Streptococcus* | *Staphylococcus aureus* | *Branhamella catarrhalis* |
| *Group G Streptococcus* | *Staphylococcus epidermidis* | *Hemophilus influenzae* |
| *Streptococcus mutans* | *Staphylococcus saprophyticus* | *Pseudomonas aeruginosa* |
| *Streptococcus intermedius* | *Corynebacterium diphtheria* | *Klebsiella pneumoniae* |
| *Streptococcus pneumoniae* | *Neisseria gonorrheae* | *Candida albicans* |
| *Streptococcus faecalis* | *Neisseria meningitidis* | *Bordetella pertussis* |

## 14. LITERATURE

1. American Academy of Pediatrics. Peter, G., éd. 1994 Red Book: Report of the Committee on Infectious Diseases. 23rd éd. Elk drove Village, IL; American Academy of Pediatrics; 1994: p. 433.
2. Bisno AL. Group A streptococcal infections and acute rheumatic fever.N. Engl. J. Med. 325: 783-793(1991).
3. Bisno AL, Pearce IA, Wall HP, Moody MD, and Stollerman GH. Contrasting epidemiology of acute rheumatic feverand acute glomerulonephritis. N. Eng. J. Med. 283: 561-565 (1970).
4. Edwards, E.A., Phillips, LA., and Suiter, W.C., Diagnosis of Group A Streptococcal Infections Directly from Throat Sécrétions, J. Clin. Micro., 15,481-483(1982).
5. Facklam RR. U.S. Dept. of Health and Human Services, PHS, CDC, Pub. No. CDC 77-13.
6. Gupta, R., Talwar, G.P. and Gupta, S.K., Rapid Antibody Capture Assay for Détection of Group A Streptococci Using Monoclonal Antibody and Colloïdal Gold Monospecific Polyvalent Antibody Conjugale, J. Immunoassay, 13,441-445(1992).
7. Heiter, BJ, Bourbeau PP, Comparison of two rapid streptococcal antigen détection assays with culture for diagnosis of streptococcal pharyngitis. J. clin. Microb. 16 : 748-53 (1995)
8. Kuttner AG and Krumwiede E. Observations on thé effect of Streptococcal upper repiratory infections on rheumatic children. a three-year study. J. Clin. Invest. 20: 273-287 (1941).
9. Lauer BA, Relier LB and Mirrell S. Effect of atmosphère and duration of incubation on primary isolation of group A Streptococci from throat cultures. J. Clin. Microb. 17:338-340 (1983)
10. Levinson, M.L. and Frank, PP., Differentiation of Group A from other Bêta Hemolytic Streptococci with Baci-tracin, J. BacterioL 69, 284-287 (1955).
11. Potter EV, Svartman M, Mohamed 1, Cox R, Poo-King T, and Earle DP. Tropical acute rheumatic fever and associated streptococcal infections compared with concurrent acute glomerulonephritis. J. Pediatr. 92: 325-333 (1978).
12. Ruoff KL, Whiley RA, Beighton D. Streptococcus. In: Manual of Clinical Microbiology, 7th éd., PR Muray, EJ Baron (eds), American Society of Microbiology, Chapter 17, pp. 283-297 (1999).
13. Ross, P.W., Throat Swabs and Swabbling Technique, The Practitioner, 207, 791-796(1971).
14. Wannamaker LW. Différences between streptococcal infections of thé throat and skin. N. Eng. J. Med. 282:78-85(1970).
15. Wannamaker LW. Changes and changing concepts in thé biology of group A Streptococci and thé epidemiology of streptococcal infections. Rev. Infect. Dis., 2: 967-973, (1979).
16. Poyart C., Bingen E., Maisonneuve P., Deschenes M., Charlier-Bret N., Boucher B. Réévaluation des Tests de Diagnostic Rapide des Angines à Streptocoque du Groupe A (Streptococcus pyogenes) (Dec. 2002). http://afssaps.santé.fr/htm/5/reactif/tdr.pdf

## 15. HISTORY

| IFU version | Amended paragraph | Details of the changes |
|---|---|---|
| V1 | All | Multilingual format and transition to V1 |
| | Content | Addition of the table of contents |
| | 15 : History | Addition of this paragraph |
| V2 | 4 : Materials provided | Remove the package insert from the content of the kit |
| | 6 : Precautions | Removal of the precautionary statement for Reagent A and addition of hazard and precautionary statements according to REACH/CLP requirements for each applicable component |

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## 14. LITTERATUR

1. American Academy of Pediatrics. Peter, G., éd. 1994 Red Book: Report of the Committee on Infectious Diseases. 23rd éd. Elk drove Village, IL; American Academy of Pediatrics; 1994: p. 433.
2. Bisno AL. Group A streptococcal infections and acute rheumatic fever.N. Engl. J. Med. 325: 783-793(1991).
3. Bisno AL, Pearce IA, Wall HP, Moody MD, and Stollerman GH. Contrasting epidemiology of acute rheumatic feverand acute glomerulonephritis. N. Eng. J. Med. 283: 561-565 (1970).
4. Edwards, E.A., Phillips, LA., and Suiter, W.C., Diagnosis of Group A Streptococcal Infections Directly from Throat Sécrétions, J. Clin. Micro., 15,481-483(1982).
5. Facklam RR. U.S. Dept. of Health and Human Services, PHS, CDC, Pub. No. CDC 77-13.
6. Gupta, R., Talwar, G.P. and Gupta, S.K., Rapid Antibody Capture Assay for Détection of Group A Streptococci Using Monoclonal Antibody and Colloïdal Gold Monospecific Polyvalent Antibody Conjugale, J. Immunoassay, 13,441-445(1992).
7. Heiter, BJ, Bourbeau PP, Comparison of two rapid streptococcal antigen détection assays with culture for diagnosis of streptococcal pharyngitis. J. clin. Microb. 16 : 748-53 (1995)
8. Kuttner AG and Krumwiede E. Observations on thé effect of Streptococcal upper repiratory infections on rheumatic children. a three-year study. J. Clin. Invest. 20: 273-287 (1941).

9.Lauer BA, Relier LB and Mirrell S. Effect of atmosphère and duration of incubation on primary isolation of group A Streptococci from throat cultures. J. Clin. Microb. 17:338-340 (1983)
10.Levinson, M.L. and Frank, PP., Differentiation of Group A from other Bêta Hemolytic Streptococci with Baci-tracin, J. BacterioL 69, 284-287 (1955).
11.Potter EV, Svartman M, Mohamed 1, Cox R, Poo-King T, and Earle DP. Tropical acute rheumatic fever and associated streptococcal infections compared with concurrent acute glomerulonephritis. J. Pediatr. 92: 325-333 (1978).
12.Ruoff KL, Whiley RA, Beighton D. Streptococcus. In: Manual of Clinical Microbiology, 7th éd., PR Muray, EJ Baron (eds), American Society of Microbiology, Chapter 17, pp. 283-297 (1999).
13.Ross, P.W., Throat Swabs and Swabbling Technique, The Practitioner, 207, 791-796(1971).
14.Wannamaker LW. Différences between streptococcal infections of thé throat and skin. N. Eng. J. Med. 282:78-85(1970).
15.Wannamaker LW. Changes and changing concepts in thé biology of group A Streptococci and thé epidemiology of streptococcal infections. Rev. Infect. Dis., 2: 967-973, (1979).
16.Poyart C., Bingen E., Maisonneuve P., Deschenes M., Charlier-Bret N., Boucher B. Réévaluation des Tests de Diagnostic Rapide des Angines à Streptocoque du Groupe A (Streptococcus pyogenes) (Dec. 2002). http://afssaps.santé.fr/htm/5/reactif/tdr.pdf.

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## 14. BIBLIOGRAFIA

1. American Academy of Pediatrics. Peter, G., éd. 1994 Red Book: Report of the Committee on Infectious Diseases. 23rd éd. Elk drove Village, IL; American Academy of Pediatrics; 1994: p. 433.
2. Bisno AL. Group A streptococcal infections and acute rheumatic fever.N. Engl. J. Med. 325: 783-793(1991).
3. Bisno AL, Pearce IA, Wall HP, Moody MD, and Stollerman GH. Contrasting epidemiology of acute rheumatic feverand acute glomerulonephritis. N. Eng. J. Med. 283: 561-565 (1970).
4. Edwards, E.A., Phillips, LA., and Suiter, W.C., Diagnosis of Group A Streptococcal Infections Directly from Throat Sécrétions, J. Clin. Micro., 15,481-483(1982).
5. Facklam RR. U.S. Dept. of Health and Human Services, PHS, CDC, Pub. No. CDC 77-13.
6. Gupta, R., Talwar, G.P. and Gupta, S.K., Rapid Antibody Capture Assay for Détection of Group A Streptococci Using Monoclonal Antibody and Colloïdal Gold Monospecific Polyvalent Antibody Conjugale, J. Immunoassay, 13,441-445(1992).
7. Heiter, BJ, Bourbeau PP, Comparison of two rapid streptococcal antigen détection assays with culture for diagnosis of streptococcal pharyngitis. J. clin. Microb. 16 : 748-53 (1995)
8. Kuttner AG and Krumwiede E. Observations on thé effect of Streptococcal upper repiratory infections on rheumatic children. a three-year study. J. Clin. Invest. 20: 273-287 (1941).
9. Lauer BA, Relier LB and Mirrell S. Effect of atmosphère and duration of incubation on primary isolation of group A Streptococci from throat cultures. J. Clin. Microb. 17:338-340 (1983)
10. Levinson, M.L. and Frank, PP., Differentiation of Group A from other Bêta Hemolytic Streptococci with Baci-tracin, J. BacterioL 69, 284-287 (1955).
11. Potter EV, Svartman M, Mohamed 1, Cox R, Poo-King T, and Earle DP. Tropical acute rheumatic fever and associated streptococcal infections compared with concurrent acute glomerulonephritis. J. Pediatr. 92: 325-333 (1978).
12. Ruoff KL, Whiley RA, Beighton D. Streptococcus. In: Manual of Clinical Microbiology, 7th éd., PR Muray, EJ Baron (eds), American Society of Microbiology, Chapter 17, pp. 283-297 (1999).
13. Ross, P.W., Throat Swabs and Swabbling Technique, The Practitioner, 207, 791-796(1971).

14. Wannamaker LW. Différences between streptococcal infections of thé throat and skin. N. Eng. J. Med. 282:78-85(1970).
15. Wannamaker LW. Changes and changing concepts in thé biology of group A Streptococci and thé epidemiology of streptococcal infections. Rev. Infect. Dis., 2: 967-973, (1979).
16. Poyart C., Bingen E., Maisonneuve P., Deschenes M., Charlier-Bret N., Boucher B. Réévaluation des Tests de Diagnostic Rapide des Angines à Streptocoque du Groupe A (Streptococcus pyogenes) (Dec. 2002). http://afssaps.santé.fr/htm/5/reactif/tdr.pdf.

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