new art laboratories

NADAL® hCG Pregnancy Test (test strips)

REF 151902N-50 / 151903N-50 / 153902N-25/ 153903N-25

de	Gebrauchsanweisung	2	da	Brugsanvisning	58
en	Instructions for use	8	no	Bruksanvisning	63
fr	Notice d'utilisation	13	nl	Gebruiksaanwijzing	68
es	Instrucciones de uso	19	ro	Instrucțiuni de utilizare	74
it	Istruzioni per l'uso	25	bg	Инструкции за употреба	79
pl	Instrukcja używania	31	el	Οδηγίες χρήσης	85
pt	Instruções de utilização	37	sk	Návod na použitie	91
cs	Návod k použití	43		Symbols	98
fi	Käyttöohjeet	48		Our International Teams	100
sv	Bruksanvisning	53			

CE IVD

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Director:
Lukas Eder

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

Version 1.1, 2025-12-12

NADAL® hCG Pregnancy Test

(Ref. 151902N-50/151903N-50/153902N-25/153903N-25)

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15. Referenzen

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update (12(6): 769-84
6. http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

Rev. 1, 2025-12-12 OM

NADAL® hCG Pregnancy Test

(Ref. 151902N-50/151903N-50/153902N-25/153903N-25)

These instructions for use contain information on the following NADAL® hCG Pregnancy Tests:

Ref.	Detection level	Packaging
151903N-50	10 mIU/mL in urine or serum (early test)	Test strips, individually packed
153902N-25	10 mIU/mL in urine or serum (early test)	Test strips in a tube
151902N-50	25 mIU/mL in urine or serum	Test strips, individually packed
153903N-25	25 mIU/mL in urine or serum	Test strips in a tube

1. Intended Use

The NADAL® hCG Pregnancy Test (151903N-50/ 153902N-25/ 151902N-50/ 153903N-25) is a lateral flow chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) at the concentration of 10 mIU/mL or 25 mIU/mL respectively in human urine or serum specimens. The test is intended for use as an aid in the early diagnosis of pregnancy (see section 12 'Limitations'). Positive results must be confirmed using (an) alternative testing method(s) (such as ultrasonography or ELISA) and clinical findings. The test procedure is not automated and requires no special training or qualification. The NADAL® hCG Pregnancy Test is designed for professional use only.

2. Introduction and Clinical Significance

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilisation and can already be detected in urine and serum in the early stage of pregnancy.

In normal pregnancies, hCG can be detected in both serum and urine as early as 7 days after conception. The hCG level continues to rise very rapidly, frequently exceeding 50 mIU/mL by the first missed menstrual period, and culminates in the 10th-12th weeks of pregnancy in the range of 30,000-290,000 mIU/mL.

The early appearance of hCG in urine and blood soon after conception, and its subsequent rapid rise in concentration, makes it an ideal marker for the early detection of pregnancy.

Elevated hCG levels comparable to the early stages of pregnancy may be associated with a hydatidiform mole or choriocarcinoma. For this reason, these possibilities should be ruled out before a pregnancy is diagnosed.


Expected values:

The reference range (serum, 97.5 percentile, beta hCG) of healthy, fertile, non-pregnant women is <3 mIU/mL. In postmenopausal women this value increases to <5 mIU/mL. However, some isolated higher values can be found in other literature. There is a sharp rise following fertilisation, so that approximately 10 days afterwards (4 days before the missed period), the values of ≥10 mIU/mL have often already been reached. The doubling time varies throughout the course of the pregnancy and is approximately 2.5 days at the beginning of the pregnancy. Due to variations in liquid absorption, there is a high individual variability for hCG concentrations in urine. Over the course of the first trimester, concentrations of 30,000 - 290,000 mIU/mL are determined, which then decrease in the 2nd and 3rd trimester to values of 45,000 - 100,000 mIU/mL. Following birth, hCG serum values decrease and are no longer detectable after approximately 11-17 days. [1], [6]

3. Test Principle

The NADAL® hCG Pregnancy Test enables the detection of human chorionic gonadotropin through the visual interpretation of colour development on the test strip. Specific anti-hCG antibodies are immobilised in the test line region (T) of the membrane. During the test, the specimen reacts with specific anti-hCG antibodies which are conjugated to coloured particles and precoated onto the conjugate pad on the test strip. The mixture then migrates along the membrane by capillary action and interacts with the reagents on the membrane. If there is a sufficient amount of hCG in the specimen, a coloured line will develop in the test line region (T) of the membrane. The presence of this coloured line indicates a positive result, whereas its absence indicates a negative result.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Materials Supplied

- 151903N-50: 50 NADAL® hCG Pregnancy test strips, individually packed
- 153902N-25: 25 NADAL® hCG Pregnancy test strips in a tube
- 151902N-50: 50 NADAL® hCG Pregnancy test strips, individually packed
- 153903N-25: 25 NADAL® hCG Pregnancy test strips in a tube
- 1 instructions for use

5. Additional Materials Required

- Specimen collection containers
- Timer

6. Storage and Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test strips are stable until the expiry date printed on the foil pouches or tubes. Test strips must remain in the sealed foil pouches or the sealed tubes until use. Reseal the tube immediately after removing a test strip. This will ensure that the remaining test strips remain dry. Once the tube has been opened, the test strips can be used for up to 90 days. Note the date when the tube was first opened and ensure that the tests

are used within this period of time. Do not freeze test kits. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

Note: All expiry dates are printed in Year-Month-Day format. 2022-06-18 indicates June 18, 2022.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the instructions for use prior to testing.
- Do not use tests beyond the expiry date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Do not use the test if the foil pouch does not contain a desiccant bag or if the desiccant bag is damaged inside the foil pouch.
- Tests are for single use only.
- Do not immerse the test strip beyond the maximum immersion line.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of in accordance with local regulations.

8. Specimen Collection, Preparation and Storage

Urine specimens:

Collect urine specimens in clean and dry plastic or glass containers which are free from any preservatives. From the first day of the missed period, hCG levels are normally so high that urine specimens collected at any time of the day may be used. However, excessive fluid intake should be avoided (dilution effect).

If the test is used at an earlier point in time, first morning urine should be used, as it generally contains the highest hCG concentration. For tests with a cut-off level of 25 mIU/mL, testing can be carried out starting from 2 days before the expected period. For tests with a cut-off level of 10 mIU/mL, the earliest point in time for testing is 4 days before the expected period. However, the risk of false negative results increases the earlier the testing is carried out (see also point 12 'Limitations').

Urine specimens exhibiting visible precipitates should be centrifuged, filtered or allowed to settle to obtain clear supernatant for testing. Only clear supernatant should be used for testing.

Do not leave specimens at room temperature for a prolonged period of time. If testing is not performed immediately, specimens can be stored at 2-8°C for up to 48 hours. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Serum specimens:

Separate serum from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum specimens can be stored at 2-8°C for up to 2 days. For long-term storage, specimens should be kept at -20°C.

Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.


If specimens are to be shipped, they should be packed in compliance with all applicable regulations for the transportation of etiologic agents.

Icteric, lipemic, haemolysed, viscous, heat-treated and contaminated specimens may lead to inaccurate test results.

9. Test Procedure

Bring tests, specimens and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test strip from the foil pouch or from the tube and use it as soon as possible. **Note:** Once the strip is removed from the tube, close it immediately afterwards. The best results will be obtained if the test is performed immediately after opening the foil pouch or tube. Label the test strip with the patient or control identification.
2. Hold the test strip vertically at the marked end and immerse it in the specimen for at least 15 seconds. In order to avoid contamination, do not touch the membrane on the test strip, and make sure not to


NADAL® hCG Pregnancy Test

(Ref. 151902N-50/151903N-50/153902N-25/153903N-25)

immerse the test strip beyond the maximum immersion line (MAX).


3. Remove the test strip from the specimen and place it on a clean and level surface.
4. Start the timer.
 As the test begins to run, you will observe a coloured liquid migrate along the membrane.
5. Wait for the coloured line(s) to appear. Read the test result after exactly 5 minutes. Do not interpret the result after more than 5 minutes.

10. Result Interpretation

Positive:

A coloured line develops in the control line region (C) and another coloured line develops in the test line region (T).

Note: The colour intensity in the test line region (T) may vary depending on the analyte concentration present in the specimen. Any shade of colour in the test line region (T) should be considered positive. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in the control line region (C). No coloured line develops in the test line region (T).

Invalid:

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Please review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test strip:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, correct procedural technique and adequate membrane wicking.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® hCG Pregnancy Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of human chorionic gonadotropin (hCG) in human urine or serum specimens only.
- Neither the quantitative value nor the rate of increase/decrease in the hCG concentration can be determined using this qualitative test.
- The accuracy of the test depends on the quality of the specimen. Inaccurate test results may occur due to improper specimen collection or storage (see section 8 'Specimen Collection and Preparation').
- Very diluted urine specimens exhibiting low specific gravity may not contain representative levels of hCG. If pregnancy is still suspected after a negative test result, a first morning urine specimen should be collected 48-72 hours later and tested.
- If tests are carried out earlier than before the first day of a missed period, there is a risk of false negative results due to individual fluctuations of hCG concentrations. According to the literature, if 10 mIU/mL tests are carried out 4 days before the expected period, only about 53-74% of pregnancies are initially detected when testing urine specimens. If pregnancy is suspected, testing should therefore be carried out again with morning urine at a later point in time. Alternatively, the hCG concentration can be determined using a more sensitive quantitative method.
- In healthy men and healthy, fertile, non-pregnant women, no hCG can be normally detected using NADAL® hCG Pregnancy Tests. Basal hCG values increase in the menopause, so that in rare cases, especially with 10 mIU/mL tests, false positive results may occur.
- Very low levels of hCG (<50 mIU/mL) are present in urine or serum shortly after the implantation of a fertilised ovum in the womb. However, since a significant number of first trimester pregnancies miscarry for natural reasons, a test result that is weakly positive should be confirmed using other clinical methods. If necessary, the test should be carried out again with the first morning urine at a later point in time. The result of the retesting should be clearly positive.
- If, after a weak positive result, the result of the retesting using the first morning urine at a later point in time is negative, a spontaneous miscarriage might be the cause. Approximately 22% of clinically undetected pregnancies end with an unnoticed miscarriage and are perceived as a slightly delayed period.
- A number of conditions other than pregnancy, including trophoblastic diseases, such as a molar pregnancy, choriocarcinoma, testicular tumors, prostate cancer and lung cancer - as well as medication containing hCG - cause elevated levels of hCG, which can lead to a positive test result.
- Positive hCG levels may be detectable for several weeks following delivery or abortion.
- Genito-urinary infections can lead to a false-positive hCG test result. An increased number of leucocytes, a basic pH, an increased amount of nitrite, erythrocytes or proteins in urine among others are indicative of genito-urinary infections. In case of suspicion, please verify this by means of urinanalysis strips or specific swab tests. Very suspicious urine specimens should be verified by laboratory testing. In the case of a positive result, no pregnancy should be diagnosed until the aforementioned possibilities have been excluded.

- An ectopic pregnancy cannot be distinguished from a normal pregnancy only on the basis of measuring hCG levels.
- The tests are calibrated for intact hCG. The proportions of intact hCG and other hCG forms in serum and urine may differ according to the pregnancy stage. Especially at the beginning of pregnancy, when hCG values are still low, this can lead, in rare cases, to deviating results for both specimen materials. In rare cases, at the beginning of pregnancy, deviating results may also be observed with tests based on the detection of beta-hCG.
- For urine specimens: In rare cases, at a more advanced stage of pregnancy (approximately week 5-8 of pregnancy), increased amounts of hCG-beta-core fragments may be present in the urine, hindering the detection reaction and leading to false negative results despite a high hCG concentration. [7]
- For serum specimens: As with all assays that utilise anti-mouse antibodies, there is a possibility of interference with anti-mouse antibodies (HAMA) in the specimen. Specimens from patients who have received preparations of monoclonal antibodies for diagnosis or therapy may contain HAMA (immunologically interfering substances). Such specimens may cause false positive or false negative test results. Heterophile antibodies in serum may also cause erroneous results.
- The test uses anti-hCG antibodies to detect hCG in urine or serum specimens. The specificity of this immune reaction rules out cross-reactivity with structurally-related hormones (hFSH, hLH and hTSH) at or above the maximum anticipated physiological levels.
- As with all diagnostic tests, all results should be interpreted by a physician in conjunction with other available clinical information.
- The NADAL® hCG Pregnancy Test only detects the presence of hCG in human urine or serum specimens and should not be used as the sole criterion for a diagnosis of pregnancy.
- If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of pregnancy.

13. Performance Characteristics

Clinical performance

Diagnostic sensitivity and specificity

The NADAL® hCG Pregnancy Tests were evaluated using clinical specimens in comparison with an ELISA method.

The results are presented in the following tables:

Ref. 151902N-50 / 153903N-25, urine/serum test (25 mIU/mL), specimen: urine

		ELISA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	210	0	210
	Negative	0	300	300
	Total	210	300	510

Diagnostic sensitivity: >99.0% (98.3% - 100.0%)*

Diagnostic specificity: >99.0% (98.8% - 100.0%)*

Overall agreement: >99.0% (99.3% - 100.0%)*

*95% confidence interval

Ref. 151902N-50 / 153903N-25, urine/serum test (25 mIU/mL), specimen: serum

		ELISA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	80	2	82
	Negative	0	259	259
	Total	80	261	341

Diagnostic sensitivity: >99.0% (98.5% - 100.0%)*

Diagnostic specificity: 99.2% (97.2% - 100.0%)*

Overall agreement: 99.4% (97.9% - 100.0%)*

*95% confidence interval

Ref. 151903N-50/ 153902N-25, urine/serum test (10 mIU/mL), specimen: urine

		ELISA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	94	0	94
	Negative	0	106	106
	Total	94	106	200

Diagnostic sensitivity: >99.0% (96.2% - 100.0%)*

Diagnostic specificity: >99.0% (96.6% - 100.0%)*

Overall agreement: >99.0% (98.2% - 100.0%)*

*95% confidence interval

Ref. 151903N-50/ 153902N-25, urine/serum test (10 mIU/mL), specimen: serum

		ELISA		
		Positive	Negative	Total
NADAL® hCG Pregnancy Test	Positive	94	0	94
	Negative	0	106	106
	Total	94	106	200

Diagnostic sensitivity: >99.0% (96.2% - 100.0%)*

Diagnostic specificity: >99.0% (96.6% - 100.0%)*

Overall agreement: >99.0% (98.2% - 100.0%)*

*95% confidence interval

Analytical performance

Detection limit

The detection limit of the NADAL® hCG Pregnancy Test is 10 mIU/mL for 151903N-50/ 153902N-25 or 25 mIU/mL for 151902N-50/ 153903N-25. The test has been standardised to the WHO International Standard.

Measuring range

Positive results were obtained when testing specimens containing hCG concentrations up to 500,000 mIU/mL. A decrease in T-line intensity was observed at concentrations >10,000 mIU/mL, indicating the onset of a prozone effect. However, hCG concentrations as high as 500,000 mIU/mL were reliably detected as positive. Thus, the measuring range of the NADAL® hCG Pregnancy Test is between 10 mIU/mL and 500,000 mIU/mL for 151903N-50/ 153902N-25 or 25 mIU/mL and 500,000 mIU/mL for 151902N-50/ 153903N-25.

Analytical specificity

Cross-reactivity study

hCG negative controls, as well as hCG positive controls (10 mIU/mL hCG for 151903N-50/ 153902N-25 and 25 mIU/mL hCG for 151902N-50/ 153903N-25) spiked with the following structurally related hormones showed no cross-reactivity when tested using the NADAL® hCG Pregnancy Tests:

hLH (300 mIU/mL), hFSH (1000 mIU/mL) and hTSH (1000 µIU/mL).

Interference study

hCG negative controls, as well as hCG positive controls spiked with the following potentially interfering substances at the concentrations listed below, showed no interference with the NADAL® hCG Pregnancy Tests:

Substance	Concentration
Paracetamol (Acetaminophen)	0.2 mg/mL
Acetylsalicylic acid	0.2 mg/mL
Albumin	20 mg/mL
Amikacin	0.2 mg/mL
Ampicillin	0.2 mg/mL
Ascorbic acid	0.2 mg/mL
Aspartame	0.2 mg/mL
Atropine	0.2 mg/mL
Atropine sulfate	0.2 mg/mL
Benzoic acid	0.2 mg/mL
Benzylpenicillin (penicillin G)	0.2 mg/mL
Bilirubin	0.01 mg/mL
Caffeine	0.1 mg/mL
Citrate	3.20%
Creatinine	0.2 mg/mL
Desoxyephedrine	0.2 mg/mL
Dextromethorphan	0.2 mg/mL
EDTA	0.2 mg/mL
Ethanol	1.0%
Gentisic acid	0.2 mg/mL
Glucose	20 mg/mL
Haemoglobin	0.01 mg/mL
Heparin	1%
Histamine	0.2 mg/mL
Methanol	1%
Methaqualone	0.2 mg/mL
Oxalic acid	0.4 mg/mL
Phendimetrazine	0.2 mg/mL
pH4-9 (positive urine specimens)	
pH4-9 (negative urine specimens)	
Quinine	0.2 mg/mL
Ranitidine	0.2 mg/mL
Sodium salicylate	0.2 mg/mL
Tetracycline	0.2 mg/mL
Tetryzoline	0.2 mg/mL
Tryptophan	0.2 mg/mL
Urea	20 mg/mL
Uric acid	0.2 mg/mL

Precision

Repeatability

Repeatability was established by testing 20 replicates of specimens at 4 concentration levels (negative, 3 mIU/mL, 10 mIU/mL and 25 mIU/mL for 151903N-50/ 153902N-25 as well as negative, 5 mIU/mL, 25 mIU/mL and 50 mIU/mL for 151902N-50/ 153903N-25) using 3 lots of the NADAL® hCG Pregnancy Tests. >99% of the specimens were correctly identified (80/80 correct tests, 95% confidence interval: 95.5% - 100%). The NADAL® hCG Pregnancy Test demonstrated acceptable repeatability.

Reproducibility

Reproducibility was established by testing 5 replicates of specimens at 4 concentration levels (negative, 3 mIU/mL, 10 mIU/mL and 25 mIU/mL for 151903N-50/ 153902N-25 as well as negative, 5 mIU/mL, 25 mIU/mL and 50 mIU/mL for 151902N-50/ 153903N-25). Testing was performed by 3 operators using 3 independent NADAL® hCG Pregnancy test lots at 3 sites on 5 separate days. >99% of the specimens were correctly identified (300/300 correct tests, 95% confidence interval: 99.8% - 100%). The NADAL® hCG Pregnancy Test demonstrated acceptable reproducibility.

14. Serious incident reporting

In the case of any serious incidents related to the performance of the NADAL® hCG Pregnancy Test, please inform nal von minden GmbH and the competent authority immediately. If still possible, do **not** dispose of the used test and the corresponding test kit components.

15. References

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update (12(6): 769-84
6. http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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15. Referencias

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff
2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136

3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update (12(6): 769-84
6. http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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15. Bibliografía

1. Ed. Lothar Thomas, Labor und Diagnose (8th Edition, 2012) TH Books, Chapter 38.3, 1859 ff

2. Butler SA et al. (2001) Detection of Early Pregnancy Forms of Human Chorionic Gonadotropin by Home Pregnancy Test Devices, Clinical Chemistry 47(12): 2131-2136
3. Lenton EA et al. (1982) Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertility and Sterility, 37(6): 773-778
4. Fritz MA et al. (1987) Doubling time of human chorionic gonadotropin (hCG) in early normal pregnancy: relationship to hCG concentration and gestational age, Fertility and Sterility, 47(4): 584-589
5. Stenman U-H et al. (2006) The classification, functions and clinical use of different isoforms of HCG, Human Reproduction Update (12(6): 769-84
6. http://americanpregnancy.org/while-pregnant/hcg-levels/from 15-8-19, 2:08 pm
7. Gronowski AM et al. (2009) False-Negative Results in Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Due to Excess hCGβ Core Fragment, Clinical Chemistry, 55 (7): 1389-1394
8. Gronowski AM et al (2009) False-Negative Results from Point-of-Care Qualitative Human Chorionic Gonadotropin (hCG) Devices Caused by Excess hCG Core Fragment Vary with Device Lot Number; Clin Chem 55:10, 1885-1886

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