new art laboratories

NADAL® Rotavirus Test (test cassette)

REF 481017

DE	Gebrauchsanweisung	2
EN	Instructions for use	5
FR	Instructions d'utilisation	8
ES	Instrucciones de uso	11
IT	Istruzioni per l'uso	14
PL	Sposób użycia	17
PT	Instruções de Utilização	20
CZ	Návod k použití	23
	Symbols	27
	Our Teams	28

IVD CE

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Moers
Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Regensburg
Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

Version 1.4, 2019-05-23

NADAL® Rotavirus Test (Ref. 481017)

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14. Referenzen

1. WILHELMI I, ROMAN E, SANCHEZ-FAUQUIER A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262
2. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38
3. Hung, T et al (1984) Waterborne outbreak of Rotavirus Diarrhoea in Adults in China caused by a Novel Rotavirus. Lancet, May 26;1(8387): 1139-1142
4. Cukor, G; Perron, DM; Hudson, R and Blacklow, NR (1984) Detection of Rotavirus in Human Stools by Using Monoclonal Antibody. J. Clin. Microl. 19: 888-892
5. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954

Rev. 4, 2019-05-23 SW

NADAL® Rotavirus Test (Ref. 481017)

1. Intended Use

The NADAL® Rotavirus Test by nal von minden is a rapid chromatographic immunoassay for the qualitative detection of rotavirus in human fecal specimens to aid in the diagnosis of rotavirus infection (rapid test). The test is designed for professional *in-vitro* diagnostic use only.

2. Introduction and Clinical Significance

Rotavirus is the most common agent responsible for acute gastroenteritis, mainly in young children.

Its discovery in 1973 and its association with infantile gastroenteritis represented a very important advancement in the study of gastroenteritis and the finding that it was not caused by an acute bacterial infection. Rotavirus is transmitted by an oral-faecal route with an incubation period of 1-3 days. Although specimen collections taken within the second and fifth day of the illness are ideal for antigen detection, the rotavirus may still be found while diarrhea continues. Rotaviral gastroenteritis may result in mortality for demographics at risk such as infants, the elderly, and immunocompromised patients.

In temperate climates, rotavirus infections occur mainly in the winter months. Endemics, as well as epidemics, affecting some thousand people have been reported. With hospitalised children suffering from acute entric disease up to 50% of the analysed specimen were positive for rotavirus. The viruses replicate in the cell nucleus and tend to be host species specifically producing a characteristic cytopathic effect (CPE). Because rotavirus is extremely difficult to culture, it is unusual to use virus isolation in diagnosing an infection. Instead, a variety of techniques have been developed to detect rotavirus in feces.

The NADAL® Rotavirus Test is a rapid chromatographic immunoassay for the qualitative detection of rotavirus in human feces specimen, providing results in 10 minutes. The test utilises antibodies specific for rotavirus to selectively detect rotavirus from human feces specimens.

3. Test Principle

The NADAL® Rotavirus Test is a qualitative, lateral flow immunoassay for the detection of rotavirus in human feces specimens.

In this assay rotavirus is detected with the aid of specific antibodies. After the addition of the sample (feces diluted with a buffer) a colour-labelled antibody specifically binds to the virus if it is present in the sample. When this complex migrates along the membrane by capillary action, it is captured with the aid of another specific antibody at the test result line region of the test. A red test result line is then generated. If no virus is present the colour labelled antibody cannot bind at the test result line region. No red test result line is formed. Therefore the presence of a coloured test result line indicates a positive result, while its absence indicates a negative result.

To serve as a procedural control, a coloured line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 individually pouched NADAL® Rotavirus test cassettes
- 10 specimen collection tubes with dilution buffer for sample collection and dilution
- 10 disposable pipettes for sample collection of extremely liquid specimen
- 1 package insert

5. Additional Materials Required

Patient:

- Stool collector. Ensure that the stool sample is collected in a way that ensures no contact with the water in the toilet bowl to avoid a dilution of the sample or contamination with detergents. On request special stool collection units can be provided by nal von minden GmbH.

Doctor's office or lab:

- Tissue paper for breaking the tip of the buffer tube
- Timer

6. Storage & Stability

Store kit as packaged either at room temperature or refrigerated (2-30°C). Under these conditions the test and the buffer are stable until the expiration date. The test must remain in the sealed pouch containing a desiccant until use. Do not freeze. Do not use beyond the expiration date.

7. Warnings and Precautions

- For *in-vitro* diagnostic use by professionals only.
- For single use. Do not reuse tests
- Do not interchange or mix reagents from different lots.
- Do not use test if its foil pouch has been damaged. Do not use tests after expiration date.
- Handle all specimens as if they contain infectious agents. Do not eat, drink or smoke in the area where the specimens or kits are handled. Protective clothing such as laboratory coats, disposable gloves and eye protection are recommended. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens in accordance with local regulations.
- The dilution buffer contains low amounts of sodium azide.
- Humidity and temperature can adversely affect results.
- The components of the test (e.g. antibodies/chemicals) do not cause any danger if the test is used according to the instructions.
- Please follow the instructions for use carefully. Please inform your patients how to collect and dilute the stool sample.

8. Specimen Collection and Preparation

Conditions for optimal sample collection

- The NADAL® Rotavirus Test is intended only for use with human fecal specimen that has been diluted in the provided buffer.
- Viral detection is improved by collecting the specimens at the onset of the symptoms. It has been reported that the maximum excretion of rotavirus in the feces of patients with gastroenteritis occurs 2-5 days after the onset of symptoms. If the specimens are collected long after the onset of diarrheic symptoms, the quantity of antigen may

not be sufficient to obtain a positive reaction or the antigens detected may not be linked to the diarrheic episode.

- Stool samples should be stored at 2-8°C immediately after collection and processed within 48 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided.

Specimen collection and preparation by the patient

In order to collect the stool sample, the patient is given one of the specimen collection tubes from the test kit and a disposable pipette. Please advise the patient to collect the stool sample in the following way:

1. Any clean and dry container / water repellent paper might be used for sample collection. Please ensure that the stool sample has no direct contact with the water of the toilet bowl to avoid a dilution or a contamination with detergents. An amount of 1-2 ml or 1-2 g stool is sufficient.

2. Transfer a small amount of the stool into the specimen collection tube:


For Solid Specimens:

Unscrew the cap of the specimen collection tube, then randomly insert the specimen collection applicator into at least 3 different sites of the faecal specimen to collect approximately 50 mg of feces (equivalent to 1/4 of a pea).


For Liquid Specimens:

If the stool is too liquid to stick to the applicator please use the provided disposable pipette. Hold the pipette vertically, aspirate some stool, and then transfer 2 drops (approximately 50 µL) into the specimen collection tube containing the dilution buffer.

3. Place the applicator back into the tube and screw the cap tightly.
4. Shake the specimen collection tube to mix the specimen and the dilution buffer. Be careful not to break the tip of the collection tube.
5. Wrap the sample in a plastic bag and store it in a cool place. Return the sample to the doctor's office within the next 24 hours.


Note:


If the patient feels uncomfortable in diluting the stool specimen into the collection tube themselves, they may also return a container with the unprocessed stool sample to the doctor's office. The transfer of the specimen into the buffer of the tube can then be carried out as described above by the personnel of the doctor's office or lab.

9. Test Procedure

1. Allow the test device and the diluted stool sample to reach room temperature (15-30°C) prior to testing.
2. Remove the test device from its pouch when ready to perform the test. The device must have room temperature to avoid condensation of moisture on the membrane. Label the device on the provided space with patient or control identification. Opened tests should be used immediately.

3. Shake the collection tube thoroughly to ensure proper mixing of the fecal sample with the extraction solution.
4. Using a piece of tissue paper, break the tip of the collection tube with a twisting motion.
5. Hold the collection tube vertically and dispense 3 drops of solution into the round sample well of the test device by applying a gentle pressure to the walls of the tube. Avoid air bubbles in the sample well or splashes of liquid into the result window.

6. Start the timer. As the test begins to work, you will see a redish coloured liquid front moving across the membrane.
7. Wait for the coloured lines to appear. The result should be read after 10 minutes. Strong positive results may be observed sooner. Do not interpret the result after more than 20 minutes.


10. Result Interpretation

For reading the test result, the colored lines that have appeared in the test result window are interpreted as follows:

Positive

Two distinct coloured lines appear on the membrane. One coloured line forms in the control line region (C) and another apparent coloured line appears in the test result line region (T).

Note:

The intensity of the colour in the test result line region (T) will vary depending on the concentration of rotavirus antigen present in the specimen. Therefore, any shade of colour in the test line region (T) indicates that Rotavirus has been detected. Do not try to quantify the amount of antigen with this qualitative assay.

Negative

Only one coloured line appears in the control line region (C). No line appears in the test line region (T). The absence of the test result line (T) indicates that no Rotavirus has been detected by the assay.

Invalid

The control line (C) fails to appear. Results from any test that has not produced a control line at the specified reading time should be discarded. Insufficient specimen volume, insufficient specimen migration, or incorrect procedural

techniques are the most likely reasons for control line failure.

Review the procedure and repeat the test with a new test device. If the presence of visible particles inhibited the migration they might be removed by centrifugation or sedimentation. Transfer a part of the sample into a tube, sediment the particles by a brief centrifugation and pipette approximately 80-120 µl of the supernatant into the sample well of a new test device. Alternatively allow the particles to settle in the upright tube and use 80-120 µl from the topside of the liquid. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Note

When fecal samples are tested, the background may appear slightly yellowish due to the colour of the fecal samples. This is acceptable as long as it does not interfere with the interpretation of test result. The test is invalid if the background fails to clear and obscures the reading of the result.

11. Quality Control

An internal procedural control is included in the test. A coloured line appearing in the control line region (C) is an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls are tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

12. Limitations

- The NADAL® Rotavirus Test (Feces) is for professional *in-vitro* diagnostic use, and should be used for the qualitative detection of rotavirus only.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of rotavirus infection.
- In order to prevent pseudo-outbreaks, positive test results in newborns must be confirmed using an alternative testing method (PCR).
- If visible blood is present in faecal specimens, false-positive results cannot be ruled out and must therefore be verified using an alternative test method.

13. Performance Characteristics

The performance of the NADAL® Rotavirus Test was evaluated in comparison with the ELISA method. A total of 242 clinical specimens collected from children and young adults were determined with both assays.

The results are summarised in the following table:

Table: NADAL® Rotavirus Test vs. ELISA

		NADAL® Rotavirus Test		
		+	-	Total
ELISA	+	79	3	82
	-	0	160	160
	Total	79	163	242

Relative Sensitivity: 96.3%
Relative Specificity: >99.9%

Overall Agreement: 98.8%

Cross-reactivity

Cross-reactivity with the following organisms has been studied at 1.0 x 10^9 organisms/ml. All the organisms were found to be negative when tested with the NADAL® Rotavirus Test.

Staphylococcus aureus	***Proteus mirabilis***	***Neisseria gonorrhoeae***
Pseudomonas aeruginosa	***Acinetobacter spp***	**Group B Streptococcus**
Enterococcus faecalis	***Salmonella choleraesius***	***Proteus vulgaris***
Group C Streptococcus	***Gardnerella vaginalis***	***Enterococcus faecium***
Klebsiella pneumoniae	***Acinetobacter calcoaceticus***	***Haemophilus influenzae***
Branhamella catarrhalis	***E. coli***	***Neisseria meningitidis***
Candida albicans	***Chlamydia trachomatis***	

14. References

1. WILHELMI I, ROMAN E, SANCHEZ-FAUQUIER A. Viruses causing gastroenteritis. Clin Microbiol Infect. April. 2003, vol.9:247-262
2. Cubitt, WD (1982) Rotavirus Infection: An Unexpected Hazard in Units Caring for the Elderly. Geriatric Medicine Today 1: 33-38
3. Hung, T et al (1984) Waterborne outbreak of Rotavirus Diarrhoea in Adults in China caused by a Novel Rotavirus. Lancet, May 26;1(8387): 1139-1142
4. Cukor, G; Perron, DM; Hudson, R and Blacklow, NR (1984) Detection of Rotavirus in Human Stools by Using Monoclonal Antibody. J. Clin. Microl. 19: 888-892
5. Tan B F et al. Pseudo-outbreak of rotavirus infection in a neonatal intensive care unit. Journal of Microbiology, Immunology and Infection 2016; 49:947e954

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